

Ubi^{select}™ Kit

Cat. No. 67-0018-001
Lot. No. 30348

Storage: cryopreserved / -20°C

FOR RESEARCH USE ONLY

NOT FOR USE IN HUMANS



PRODUCT DESCRIPTION Page 1 of 5

Ubi^{select} Kit Utility

The Ubi^{select} Kit is designed to facilitate the isolation and identification of your ubiquitylated proteins of interest from a transgenic mouse-derived embryonic fibroblast (MEF) cell line expressing ubiquitin tagged with an N-terminal biotinylation signal (which is biotinylated by co-expression of the BirA enzyme). By employing a powerful capture system which utilises the transgenic MEF cell line expressing ubiquitin - which is biotinylated in the cell - and the high affinity biotin-avidin interaction alongside a high stringency washing protocol you can capture your ubiquitylated target of interest without any associated ubiquitin binding proteins and/or ubiquitylated protein binding proteins. The Ubi^{select} Kit contains an aliquot of each of the control Ubi^{select}-MEF_BirA (P0) and the Ubi^{select}-MEF_Bio-Ub (P0) cell lines and an anti-ubiquitin conjugate antibody for the detection of ubiquitylated proteins. Eluted captured products may be analysed by Western blotting using the ubiquitin conjugate specific antibody, an anti-biotin antibody or an antibody specific to your protein of interest.

Background

The post-translational modification of proteins by ubiquitin is involved in a wide range of cellular processes (Kirkin and Dikic, 2007). In any given cell the proportion of protein modified by ubiquitin is very small therefore it has been difficult to isolate and identify this post translational modification from mammalian whole cell lysates. Ubiq-

Kit Components

Product	Cat. No.	Lot No.	Amount	Storage
Ubi ^{select} ™-MEF_BirA (P0)	66-5010-001	30346	2.5x10 ⁶ cells/ml	cryopreserved
Ubi ^{select} ™-MEF_Bio-Ub (P0)	66-5011-001	30347	2.5x10 ⁶ cells/ml	cryopreserved
Mono and Polyubiquitylated conjugates, mAb (FK2) HRP linked	68-0122-025	30125	25 µg	-20°C

Cell Culture Characteristics

Species: mouse

Passage Number: passage 0

Source: mouse embryonic fibroblasts

Storage: cryopreserved

Culture Medium: DMEM (HG), 200mM L-glutamine, 100mM sodium pyruvate, 10% fetal bovine serum (FBS), biotin 0.2 µg/mL, 1x pen/strep

Quality Assurance:
Mycoplasma: not detected
Morphology: fibroblast
Viability (%): ≥95

Growth Mode: adherent

Ubi^{select} Cell Culture Protocol

Cell culture procedure for Ubi^{select}-MEF_BirA (P0) and Ubi^{select}-MEF_Bio-Ub (P0)

Thaw the cells by swirling the vial in a 37°C waterbath. Place 5 ml of media into a T25 flask and pipette cells into the flask containing the media. Place the T25 flask containing the cells into a 37°C 5% CO₂ incubator overnight and re-feed the next morning with fresh media.

When BirA and Bio-Ub MEF cells reach confluency, passage as follows.

1. Aspirate culture media from the T25 flasks and rinse cells with 5 ml PBS.
2. Aspirate PBS from the T25 flask and add 1 ml 0.05% Trypsin-EDTA (1x). Incubate the cells in a 37°C 5% CO₂ incubator for a few minutes until the MEF cells begin to detach then knock the side of the flask to fully detach the cells.
3. Add 4 ml of media, collect in a 15 ml tube and count a sample of the cells. Seed MEF cells in T75 flasks at 1.2 x 10⁴ cells/cm² (i.e. 9x10⁵ cells per T75 flask) then place in a 37°C 5% CO₂ incubator and grow until confluency (P1).
4. Re-feed MEF cells every 3-4 days where necessary. When the BirA and BioUb MEFs have reached confluency, passage for a second time (P2) as above into T75 flasks, seeding at 1.2 x 10⁴ cells/cm² (i.e. 9 x 10⁵ cells per T75 flask).



www.ubiquigent.com
Dundee, Scotland, UK

ORDERS / SALES SUPPORT

International: +1-617-245-0020
US Toll-Free: 1-888-4E1E2E3 (1-888-431-3233)
Email: sales.support@ubiquigent.com

UK HQ and TECHNICAL SUPPORT

International: +44 (0) 1382 381147 (9AM-5PM UTC)
US/Canada: +1-617-245-0020 (9AM-5PM UTC)
Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2014. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0

Cat. No. 67-0018-001
Lot. No. 30348

Storage: cryopreserved / -20°C

FOR RESEARCH USE ONLY

NOT FOR USE IN HUMANS



PRODUCT DESCRIPTION Page 2 of 5

Ubiquitin proteomics remains challenging even though the sensitivity of Mass Spectrometry (MS) has improved dramatically through the use of innovative techniques (Sylvestersen *et al.*, 2013). Various approaches employing tagged ubiquitin or ubiquitin-like molecules have been used with varying degrees of success (Peng *et al.*, 2003; Tirard *et al.*, 2012; Tsirigotis *et al.*, 2001). In recent years, the isolation of ubiquitylated proteins from neurons of *Drosophila melanogaster* using a tagged ubiquitin with a 15 amino acid long biotin-accepting peptide has been described (Franco *et al.*, 2011). This was the first demonstration that proteomics could be used to identify neuronal targets of the ubiquitin-proteasome system. This novel technique allowed for the isolation and enrichment of ubiquitin conjugates from neurons using a relatively small sample up to levels that allowed direct detection by MS and Western blotting. In addition, di-glycine signatures indicative of the ubiquitin attachment sites could also be detected on ubiquitin conjugates. Where antibodies were available for specific UPS substrates it was also possible to determine whether a substrate was mono- or polyubiquitylated (Franco *et al.*, 2011). Using a similar approach, hemizygous Biotinylated-Ubiquitin (Bio-Ub) and control BirA transgenic mouse models have been created which express either three moieties of biotin-accepting ubiquitin plus the *E. coli* enzyme BirA or the *E. coli* enzyme BirA enzyme alone (Lectez *et al.*, 2014). The

Continued on page 3

Ubiselect Cell Culture and Capture Protocol

Continued from page 2

- At passage 3 (P3), plate BirA and Bio-Ub MEFs onto 100 mm dishes at a seeding density of 6×10^5 cells per dish. Harvest BirA and Bio-Ub MEFs when they reach confluency or as necessary.

Ubiselect Capture Protocol*

- To harvest BirA and Bio-Ub MEFs, remove media from the 100 mm BirA MEF and Bio-Ub MEF cell dishes, wash with 5 ml ice cold PBS and aspirate.
- Collect the cells from each 100 mm dish by adding 200 μ l of Cell Lysis Buffer to each dish, using a cell scraper as necessary.
- Prepare separate MEF cell lysates from P3 BirA and Bio-Ub MEF cells. Homogenise using a needle and syringe and clarify the lysate by centrifuging at 4000 rpm for 10 minutes. Collect the supernatant and determine the protein concentration using standard methods (Input Sample). In the protocol below, a final volume of 2.5 ml of Input Sample derived from each cell line has been assumed.
- Keep aside Input Sample for both BirA and Bio-Ub MEF cell lysates to run on your SDS-Page gels for Western blotting.
- Equilibrate two PD10 columns by allowing 25 ml of Binding Buffer (-NEM) to flow through.
- Add 2.5 ml of P3 BirA and Bio-Ub MEF Input Samples separately to the two equilibrated PD10 columns.
- Add 3.5 ml of Binding Buffer (+NEM) to the two PD10 columns and separately collect the eluate for the BirA and Bio-Ub MEF cell lysates into 15 ml falcon tubes containing 250 μ l of 25x protease inhibitor (PI) mixture (prepare 25x PI by adding one Protease Inhibitor tablet (Roche Cat# 11836170001) to 0.4ml of Binding Buffer (-NEM)) and 250 μ l of NeutrAvidin agarose bead slurry. To prepare the NeutrAvidin bead slurry, wash the required amount of beads three times in Binding Buffer (-NEM) then generate a 50:50 bead:Binding Buffer (-NEM) slurry.
- Incubate the 4 ml of eluate plus PI mix and NeutrAvidin agarose bead suspension on a rolling platform for 40 minutes at RT then a further 2 hrs and 20 mins at 4°C.
- Wash beads in 10 ml WB1 x2 washes, centrifuge for 1 minute at 4000 rpm, discarding the supernatant after each wash.

Continued on page 3



www.ubiquigent.com
Dundee, Scotland, UK

ORDERS / SALES SUPPORT

International: +1-617-245-0020
US Toll-Free: 1-888-4E1E2E3 (1-888-431-3233)
Email: sales.support@ubiquigent.com

UK HQ and TECHNICAL SUPPORT

International: +44 (0) 1382 381147 (9AM-5PM UTC)
US/Canada: +1-617-245-0020 (9AM-5PM UTC)
Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2014. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0

Cat. No. 67-0018-001
Lot. No. 30348

Storage: cryopreserved / -20°C

FOR RESEARCH USE ONLY

NOT FOR USE IN HUMANS



PRODUCT DESCRIPTION Page 3 of 5

Continued from page 2

Ubiselect kit includes an aliquot of both the Ubiselect-MEF_Bio-Ub (P0) (Cat# 66-5011-001) and the control Ubiselect-MEF_BirA (P0) (Cat# 66-5010-001) cell lines and an anti-ubiquitin conjugate antibody for the detection of ubiquitylated proteins; see Application Note 004 for details. Ubiselect-MEF_Bio-Ub (P0) and Ubiselect-MEF_BirA (P0) have been derived from 13.5 day old embryos from these transgenic mouse lines by homogenisation and trypsinisation of the embryos minus the head and liver (Lectez *et al.*, 2014). By using cell and tissue lysates derived from Bio-Ub transgenic mice, it has been possible to isolate ubiquitylated proteins only (Lectez *et al.*, 2014). The specific isolation of ubiquitylated proteins is achieved through the high affinity interaction between biotin and avidin and high stringency washing which allows for the removal of interacting partners of ubiquitylated proteins which would otherwise predominate under the milder washing conditions of alternative ubiquitin capture systems. By utilising Bio-Ub expressing liver lysates, the isolation of ubiquitin conjugated proteins (both thioester and isopeptide linked) has been demonstrated, therefore this is the first mammalian system for ubiquitin proteomics which allows for the direct detection of ubiquitin conjugates by MS and Western blotting. Isolation of ubiquitylated proteins from Bio-Ub Mouse Embryonic Fibroblast (MEF) cell lines has also been demonstrated (See Application Note 004)

Continued on page 4

Ubiselect Cell Culture and Capture Protocol

Continued from page 2

Continue with the following wash steps centrifuging and carefully discarding the supernatant as above in-between each wash step.

- 10 ml WB2 x3 washes
- 10 ml WB3 x1 wash
- 10 ml WB4 x3 washes
- 10 ml WB1 x1 wash
- 10 ml WB5 x1 wash
- 10 ml WB6 x3 washes
16. Add 100 µl of Elution Buffer to each lot of beads, transfer to a 1.5 ml tube and incubate for 5 min at 95 °C.
17. Add the sample to a Vivaclear Mini 0.8-µm PES microcentrifuge filter unit and centrifuge for 1 minute at 6000 rpm. The recovered sample is the Output Sample.
18. The Output Sample (captured biotinylated ubiquitin - both in its free form and conjugated to ubiquitylated proteins) is now ready for analysis. For analysis by Western Blotting you may wish to run 1/20 of the Output Sample per lane by SDS-PAGE and an aliquot of the Input Sample alongside (1/20 or an amount to be determined empirically) for comparison (pre and post-purification). Analyse both your BirA and Bio-Ub MEF Input and Output Samples for comparisons (see Application Note 004 for an example).
19. For the detection of biotinylated ubiquitin conjugates by Western blotting use the Mono and Polyubiquitylated conjugates, mAb (FK2) HRP linked antibody supplied with the Ubiselect kit. The incubation conditions for this antibody are as follows: Block the Western blot membrane in 5% Milk/TBST for 1 hour at room temperature then rinse the membrane in TSBT to remove any residual milk. Then incubate the Western blot membrane with Mono and Polyubiquitylated conjugates, mAb (FK2) HRP linked antibody (1:4000) in 1% BSA/TBST for 2 hours at room temperature, wash for 5 minutes in TBST, repeat three times and develop the Western Blot membrane accordingly (for example by ECL).

*NB: This protocol has been optimised to enable the selective capture and purification of ubiquitylated proteins *only* away from ubiquitylated protein binding proteins and ubiquitin and ubiquitin chain binding proteins. If you wish to also capture such binding proteins then adjust the formulation of the Cell Lysis, Binding Buffers and Wash Buffers (WB) and the protocol according to your stringency requirements (eg through the adjustment in the binding and wash step configurations including the concentration or removal of SDS, guanidine and/or urea).

Continued on page 4



www.ubiquigent.com
Dundee, Scotland, UK

ORDERS / SALES SUPPORT

International: +1-617-245-0003
US Toll-Free: 1-888-4E1E2E3 (1-888-431-3233)
Email: sales.support@ubiquigent.com

UK HQ and TECHNICAL SUPPORT

International: +44 (0) 1382 381147 (9AM-5PM UTC)
US/Canada: +1-617-245-0020 (9AM-5PM UTC)
Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2014. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0

Cat. No. 67-0018-001
Lot. No. 30348

Storage: cryopreserved / -20°C

FOR RESEARCH USE ONLY

NOT FOR USE IN HUMANS

PRODUCT DESCRIPTION Page 4 of 5

Continued from page 3

and it is anticipated that this technology will enable researchers to gain a better understanding of ubiquitin system signalling.

References:

Franco M, (2011) A novel strategy to isolate ubiquitin conjugates reveals wide role for ubiquitination during neural development. *Mol Cell Proteomics* **10**, M110 002188.

Kirkin V and Dikic I. (2007) Role of ubiquitin- and Ubl-binding proteins in cell signaling. *Curr Opin Cell Biol* **19**, 199-205.

Lectez B, Migotti R, Lee SY, Ramirez J, Beraza N, Mansfield B, Sutherland JD, Martinez-Chantar ML, Dittmar G, Mayor U. (2014) Ubiquitin Profiling in Liver Using a Transgenic Mouse with Biotinylated Ubiquitin. *J Proteome Res* **13**, 3016-3026.


Peng J, Schwartz D, Elias J E, Thoreen C C, Cheng D, Marsischky G, Roelofs J, Finley D, Gygi SP (2003) A proteomics approach to understanding protein ubiquitination. *Nat Biotechnol* **21**, 921-926.

Sylvestersen K B, Young C, Nielsen M L. (2013) Advances in characterizing ubiquitylation sites by mass spectrometry. *Curr Opin Chem Biol* **17**, 49-58.

Tirard M, Hsiao H H, Nikolov M, Urlaub H, Melchior F, Brose N. (2012) In vivo localization and identification of SUMOylated proteins in the brain of His6-HA-SUMO1 knock-in mice. *Proc Natl Acad Sci U S A* **109**, 21122-21127.

Tsirigotis M, Thurig S, Dube M, Vanderhyden BC, Zhang M, Gray DA (2001) Analysis of ubiquitination in vivo using a transgenic mouse model *Biotechniques* **31**, 120-126.

Capture Protocol Flow Diagram



www.ubiquigent.com
 Dundee, Scotland, UK

ORDERS / SALES SUPPORT

International: +1-617-245-0003
US Toll-Free: 1-888-4E1E2E3 (1-888-431-3233)
Email: sales.support@ubiquigent.com

UK HQ and TECHNICAL SUPPORT

International: +44 (0) 1382 381147 (9AM-5PM UTC)
US/Canada: +1-617-245-0020 (9AM-5PM UTC)
Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2014. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0



Cat. No. 67-0018-001
Lot. No. 30348

Storage: cryopreserved / -20°C

FOR RESEARCH USE ONLY

NOT FOR USE IN HUMANS

PRODUCT DESCRIPTION Page 5 of 5

Buffers and Reagents

Binding Buffer (+NEM)

3M urea, 1M NaCl, 0.25% SDS and 50 mM N-ethylmaleimide in PBS

Binding Buffer (-NEM)

3M urea, 1M NaCl and 0.25% SDS in PBS.

Cell Lysis Buffer

8M urea, 1% SDS and 50mM N-ethylmaleimide in PBS, (add one protease inhibitor tablet (Roche Cat# 11836170001) to 10ml of Cell Lysis Buffer).

Elution Buffer

4× Laemmli buffer including 100 mM DTT for reducing Elution Buffer and minus DTT for non-reducing Elution Buffer

TBS Tween (TBST)

50mM Tris (pH7.6), 150mM NaCl, 0.05% Tween 20

WB1

8M urea and 0.25% SDS in PBS

WB2

6M guanidine HCl in PBS

WB3

6.4M urea, 1M NaCl and 0.2% SDS in PBS

WB4

4M urea, 1M NaCl, 10% isopropanol, 10% ethanol and 0.2% SDS in PBS

WB5

8M urea and 1% SDS in PBS

WB6

2% SDS in PBS

Reagents not included in this kit:

- High capacity NeutrAvidin agarose (Thermo Scientific Cat# 29202)
- Anti-BirA antibody (Sigma Cat# GW20013F)
- Anti-biotin (D5A7) rabbit mAb; HRP conjugate (Cell Signalling Cat# 5571)
- Urea (Sigma, Cat# U0631)
- Sodium Dodecyl Sulfate (SDS) (Sigma, Cat# L6026)
- Guanidine HCl (Sigma, Cat# G3272)
- N-Ethylmaleimide (NEM) (Thermo Fisher, Cat# 23030)
- Protease inhibitor cocktail-cOmplete mini-EDTA free (Roche, Cat# 11836170001)
- Phosphate Buffered Saline (Sigma, Cat# P4417)
- PD-10 Columns (GE Healthcare Cat# 17-0851-01)
- Vivaclear Mini 0.8-µm PES microcentrifuge

Limited Use Licence

(in respect of product Cat# 67-0018-001)

The technology used to derive the cell line(s) described in this document is owned by CIC bioGUNE and exclusively licenced for commercialization by Ubiquigent Ltd.

The Biotin AviTag technology is covered by U.S. Pat. Nos: 5,723,584; 5,874,239; 5,932,433; and 6,265,552 and includes any and all materials, methods, kits and related derivatives claimed by these patents. The purchase of the UBIQUIGENT's AviTagged products confers to the purchaser the limited right to use the AviTag™ technology for research use.

Commercial use of the AviTag™ technology to manufacture a commercial product, or use of the AviTag™ technology to facilitate or advance research which will be applied to the development of a commercial product requires a license from Avidity, LLC. Examples of Commercial use include, but are not limited to biosensors, diagnostics, therapeutic applications, proximity assays, and drug screening assays. Users desiring a commercial license should contact Avidity, LLC in the US @ 720-859-6111 or by email to Lansing@avidity.com.



www.ubiquigent.com
Dundee, Scotland, UK

ORDERS / SALES SUPPORT

International: +1-617-245-0003
US Toll-Free: 1-888-4E1E2E3 (1-888-431-3233)
Email: sales.support@ubiquigent.com

UK HQ and TECHNICAL SUPPORT

International: +44 (0) 1382 381147 (9AM-5PM UTC)
US/Canada: +1-617-245-0020 (9AM-5PM UTC)
Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2014. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0