For screening Enterohemorrhagic *Escherichia coli* O157:H7/-
CLIG Agar (cellobiose added LIG agar medium)

Characteristic:

CLIG Agar enables estimation and confirmation of *E. coli* by confirmation of acid production from lactose and cellobiose, indole and β -glucuronidase production, and is useful for screening for enterohemorrhagic *Escherichia coli* O157:H7/-.

Composition (1 L = 41.0 g)

Casein peptone	7.5 g	MUG	0.02 g
Meat peptone	2.5 g	Sodium chloride	5.0 g
Lactose	1.0 g	Phenol red	0.025 g
Cellobiose	10.0 g	Agar powder	14.9 g
Tryptophan	0.1 g		
Purified water	1,000 ml	pH 7.4 \pm	

Preparation method:

41.0 g of this product is added to 1,000 mL of purified water, and dissolved by heating. 2.5 to 3.0 mL of this agar medium is poured into each small size of test tube, and high-pressure sterilization is performed at 115 degrees for 15 minutes. After sterilization, solidify it as semi-slant medium such as TSI agar medium etc.

Test method:

Inoculate the test bacterial strains in the bottom part and on the slant part, and incubate at 35 to 37 degrees for 18 to 24 hours.

Caution: Do not keep it with a tight stopper (cap) during culture.

Interpretation of results:

	Enterohemorrhagic <i>E. coli</i> O157:H7/-	Other <i>E. coli</i>
Bottom part	Change to yellow	Change to yellow
Slant part	Change to red	Change to red
Cytochrome oxidase	Negative	Negative
Indole production *1)	Positive	Positive
Blue fluorescence under UV irradiation *2)	Negative	Positive

*1): Testing by indole spot test, *2): Ultraviolet at 365 nm

After confirming the above properties, conduct agglutination test using anti-O157 serum.

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COSMO BIO Co., LTD.