

For research use only. Not for clinical diagnosis.

Catalog No. BAM-02-042-EX

E. coli Single-stranded DNA Binding Protein (SSB)

BACKGROUND

E. coli single-stranded DNA binding protein (SSB) binds to single-stranded DNA with high specificity (1, 2). It is involved in DNA replication and recombination *in vivo*. The SSB gene was expressed as the recombinant protein in *E.coli* highly purified. The molecular mass is 18.9 kDa.

Applications: 1) Functional single-stranded DNA-binding protein for studying DNA replication and

recombination

2) Enhancement of the specificity and yield of PCR

Size: 200 μg

Form: 5 mg/ml in 20mM Tris-HCl (pH 7.6), 200mM NaCl, 1mM dithiothreitol, 1mM EDTA, 50%

glycerol

Quality: Greater than 95% purity as determined by SDS-PAGE (CBB staining)

The absence of endonucleases and exonucleases was confirmed.

Data Link: Swiss-Prot <u>P0AGE0</u>

Storage: -20°C

References: 1) Krauss,G. et al. (1981) "Escherichia coli single-strand deoxyribonucleic acid binding protein: Stability,

specificity, and kinetics of complexes with oligonucleotides and deoxyribonucleic acid." Biochemistry 20:

5346-5352 PMID: <u>7028102</u>

2) Weiner, J.H. et al. (1975) "The deoxyribonucleic acid unwinding protein of Escherichia coli. Properties and

functions in replication." J.Biol.Chem. **250**:1972-1980 PMID: $\underline{1090613}$



E. coli Single-stranded DNA Binding Protein (SSB)




Fig.1 SDS-PAGE of E.coli SSB protein




Fig.2 Binding activity to single-stranded DNA 0.02 ug/ul of M13mp18ssDNA was incubated with 0 (lane0), 0.025 (lane1), 0.05 (lane2), and 0.1 (lane3) ug/ul of SSB at 37°C for 30 min. and then 10ul aliquot was subjected to electrophoresis in agarose.

Related Products

BAM-02-042-5-EX

E. coli Single-stranded DNA Binding Protein (SSB) 1mg (5 x 200 μg)

For research use only. Not for clinical diagnosis.

Manufactured by BioAcademia,Inc.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp e-mail: export@cosmobio.co.jp

[Outside Japan] Phone: +81-3-5632-9617 [国内連絡先] Phone: +81-3-5632-9610 FAX: +81-3-5632-9618 FAX: +81-3-5632-9619