

Taq Premix with Dyes

02-102 100 reactions 02-102-5 500 reactions

Description

Taq Premix with Dye is an optimized ready-to-use solution containing *Taq* DNA Polymerase, dNTPs, MgCl₂, KCl and stabilizers. It is ideally suited to routine PCR applications from templates including

pure DNA solutions bacterial colonies and cDNA products. | Composition of PCR reaction Mixture (total 50µl)

	=	•
Applications	Taq Premix with Dye	25µl
PCR	Template	<500ng
Primer Extension	Primer 1	$0.2\sim1.0\mu M$ (final conc.)
·Colony PCR	Primer 2	$0.2\sim1.0\mu M$ (final conc.)
	Sterile distilled wate	r up to 50μl

Concentration: One unit is defined as the amount of enzyme that can incorporate 10nmols of total dNTPs into an acid-insoluble material in 30 minutes at 74°C when activated salmon sperm DNA was used as template / primer.

Quality Assurance: Greater than 95% purity as determined by SDS-PAGE (CBBstaining). The absence of endonuclease and exonucleases was confirmed.

PCR product: PCR products have one A added at the 3'-terminus. Thus, the PCR product can be used directly for cloning into a T-vector. Additionally, it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

PCR Test: Good amplification result was obtained in PCR reaction using λDNA as a template (Fig.1). Premix composition: 10 mM Tris-HCl, 50 mM MgCl₂, 0.2 mM dNTPs, 5 % Glycerol, 0.08 % NP-40, 0.05 % Tween-20, 25 units/ml Taq DNA Polymerase, pH 8.6 @ 25°C, 0.0005 % BPB, 0.001% XC

Storage Temperature: - 20°C M 1 2 3 4

Fig. 1 Amplification of λ DNA

Lane M: marker

PCR条件

98° C 10sec

57° C 30sec 25cycles

72° C 8min.

(2min in the case of 2kb DNA.)

Notes: Repeated freezing and thawing may decrease enzyme activity. Once thawed, aliquot into PCR tubes and store at -20°C. (For 50 µlPCR reactions, dispense 25 µlinto each tube.)

If you store this product at 4°C, please use it within 3 months.



TOKYO JAPAN

TEL: (81)3-5632-9617

FAX: (81)3-5632-9618
e-mail: export@cosmobio.co.jp

Fig.2