

### POLYCLONAL ANTIBODY

### Catalog No. TIP-SN-P07

# Anti α-Synuclein (61-70)

### BACKGROUND

Alpha-Synuclein, a 140-amino acid protein abundantly expressed in presynaptic terminals, is known as a component of intraneuronal or glial inclusions observed in cases of Parkinson's disease (PD), Dementia with Lewy bodies (DLB) and Multiple system atrophy (MSA). Alpha-Synuclein is a natively unfolded protein, however, fibrillization or conformational change of Alpha-synuclein is central in the pathogenesis of Alpha-synucleinopathies. The amino-terminal region of Alpha-synuclein consists of seven imperfect repeats, each 11 amino acids in length, with the consensus sequence KTKEGV. The repeats partially overlap with a hydrophobic region (amino acids 61-95). The carboxy-terminal region (amino acids 96-140) is negatively charged. These antibodies are powerful tools for biochemical and IHC analyses of neurodegenerative diseases and for evaluation of conformational changes of Alpha-synuclein.

Product type	Primary antibody			
Immunogen	EQVTNVGGAVC (α-syn 61-70)			
Raised in	Rabbit (New Zealand White)			
Source	Anti-serum			
Purification	-			
Form	Liquid. Anti-serum			
Concentration	-			
Volume	50 uL			
Label	Unlabeled			
Specificity	α-Syn			
Cross reactivity	Human, Mouse			
Storage	Store below -20°C. (below -70°C for prolonged storage).			
	Aliquot to avoid cycles of freeze/thaw.			

Application notes Recommended dilutions	•ELISA : 1/500-1/2000 •Western Blot : 1/500-1/1000 •Immunohistochemistry :1/500-1/2000
	Other applications have not been tested. Optimal dilutions/ concentrations should be determined by the end user.
References	1)Masami Masuda, et al. Inhibition of a-synuclein fibril assembly by small molecules: Analysis using epitope-specific antibodies. FEBS Letters (2009) 583, 787-791. PMID 19183551



## **RELATED PRODUCTS:**

Product Name	Quantity	Maker	Cat#
Anti α-Synuclein (1-10)	50 uL	CAC	TIP-SN-P01
Anti α-Synuclein (11-20)	50 uL	CAC	TIP-SN-P02
Anti α-Synuclein (21-30)	50 uL	CAC	TIP-SN-P03
Anti α-Synuclein (31-40)	50 uL	CAC	TIP-SN-P04
Anti α-Synuclein (41-50)	50 uL	CAC	TIP-SN-P05
Anti α-Synuclein (51-60)	50 uL	CAC	TIP-SN-P06
Anti α-Synuclein (75-91)	50 uL	CAC	TIP-SN-P08
Anti α-Synuclein (131-140)	50 uL	CAC	TIP-SN-P09
Anti α-Synuclein (9 Antibodies Set)	9*10 uL	CAC	TIP-SN-SET

For research use only, Not for diagnostic use.

