



*For research use only. Not for clinical diagnosis.*

**Catalog No. CK026**

## **Sperm diluent for in vitro fertilization (SEM-5X)**

### **BACKGROUND**

This can be used for the washing and dilution of sperm used in insemination and the preparation of isotonic percoll solution. This is sterilized ready to use.

<b>Storage</b>	Stored at 4°C in a dark room. Do not freeze.
<b>Package</b>	25 ml
<b>Expiration date</b>	The expiration date is shown on the label of the medium bottle (6 months after the medium preparation)
<b>Note</b>	SEM-5X is for research use only. Do not use in animals, humans, or diagnostic procedures.

*For research use only. Not for clinical diagnosis.*



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619





For research use only. Not for clinical diagnosis.

Catalog No. CK026

媒精用精液希釈液 (5 倍液)  
(SEM-5X)

媒精時の精液の洗浄、希釈に用いる等張パーコール液の調製に使用できます。滅菌されておりますのでそのままご利用になれます。

保存法	本回収液は、冷暗所 (2~5°C) に保存して下さい。凍結厳禁。
包装	25 ml
有効期間	製造後 6 ヶ月



人と科学のステキな未来へ

コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル  
URL : <http://www.cosmobio.co.jp/>

● 営業部 (お問い合わせ)  
TEL : (03) 5632-9610 FAX : (03) 5632-9619  
TEL : (03) 5632-9620



コスモ・バイオ株式会社