

Fluoresceinamine-labeled Sodium Chondroitin Sulfate E (Squid Cartilage)

Product code: CSR-FACS-E1

Volume: 3 mL (1 mg/mL, in phosphate buffered saline (PBS) (-))

Appearance: yellow green solution

Source of sodium chondroitin sulfate E: Squid Cartilage CAS number of sodium chondroitin sulfate: 12678-07-8

Fluorescent probe: Fluoresceinamine

CAS number of fluorescent probe: 3326-34-9

Outline: Chondroitin sulfate (CS) is a sulfated glycosaminoglycan composed of repeating disaccharide units of N-acetyl-D-galactosamine (GalNAc) and D-glucuronic acid (GlcUA). CS is abundant in cartilage and exists as unbranched polysaccharide chains covalently linked to the protein core of proteoglycans. Chondroitin sulfate E (CS-E) contains approximately 60 % of disaccharide units with two sulfate groups as 4-O- and 6-O-sulfation of GalNAc (E structure unit). This product is prepared by the fluorescent labeling of CS-E derived from squid cartilage according to the method of Ogamo et al. 1). Fluoresceinamine molecules are chemically attached to carboxyl groups of the GlucUA of CS. The endotoxin content is in accordance with the product specifications. The excitation wavelength is 490 ~ 500 nm and the emission wavelength is 515~525 nm. The enclosed Certification of Analysis lists actual values for product specifications.

Handling precautions:

- 1) Protect from light as much as possible. Product can be used at room temperature when protected from strong light.
- 2) After thawing, gently agitate the vial before use.
- 3) Store protected from light at -20°C or below. We recommend storing in aliquots appropriate for anticipated usage.
- 4) Fluorescence intensity varies with pH of the solution and is lower under acidic conditions. Note the pH of the sample solution when measuring fluorescence intensity.
- 5) This product is for research use only. Not for use in diagnostic or therapeutic procedures.

Reference:

1) Ogamo A et al.: Carbohydr. Res., **105**, 69 (1982)

For research use only. Not for clinical diagnosis.

Manufactured by





Cosmo Bio Co., Ltd.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp e-mail: export@cosmobio.co.jp

[Outside Japan] Phone: +81-3-5632-9617

[国内連絡先] Phone: +81-3-5632-9610 FAX: +81-3-5632-9618 FAX: +81-3-5632-9619