



Fluoresceinamine -labeled Sodium Hyaluronate (L1)

(The average molecular weight: $100 \times 10^3 \sim 300 \times 10^3$ Da)

Product code: CSR-FAHA-L1

Volume: 3 mL (1 mg/mL, in phosphate buffered saline (PBS) (-))

Appearance: yellow green solution

Source of sodium hyaluronate: *Streptococcus* sp.

CAS number of sodium hyaluronate: 9067-32-7

Fluorescent probe: Fluoresceinamine

CAS number of fluorescent probe: 3326-34-9


Outline: Hyaluronan (HA) is a glycosaminoglycan composed of repeating disaccharide units of N-acetyl-D-glucosamine (GlcNAc) and D-glucuronic acid (GlcUA). HA is abundant in synovial fluid, skin, umbilical cord, and vitreous body exists as unbranched polysaccharide chains. This product is prepared by the fluorescent labeling of HA according to the method of Ogamo et al.¹⁾. Fluoresceinamine molecules are chemically attached to carboxyl groups of the GlcUA of HA. This solution is dissolved in PBS (-) and sterilized by filtration. The endotoxin content is in accordance with the product specifications. The excitation wavelength is 490~500 nm and the emission wavelength is 515~525 nm. The enclosed Certification of Analysis lists actual values for product specifications.

Handling precautions:

- 1) Although degradation of FAHA-L1 by light exposure is not so much as FAHA-H1 or FAHA-M1, protect from light as much as possible. Product can be used at room temperature when protected from strong light.
- 2) After thawing, gently agitate the vial before use.
- 3) Store protected from light at -20°C or below. We recommend storing in aliquots appropriate for anticipated usage.
- 4) Fluorescence intensity varies with pH of the solution and is lower under acidic conditions. Note the pH of the sample solution when measuring fluorescence intensity.
- 5) This product is for research use only. Not for use in diagnostic or therapeutic procedures.

Reference:

- 1) Ogamo A et al.: *Carbohydr. Res.*, **105**, 69 (1982)

Manufactured by  PG
Research



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

Phone : +81-3-5632-9617 FAX : +81-3-5632-9618