




# 11-Oxo-testosterone (11-KT)

**CAS Number** : 564-35-2

**Formula** :  $C_{19}H_{26}O_3$



**Molecular Weight** : 302.41

**Melting point** : 184 - 187 °C

**Purity** : Not less than 98.0 %

**Appearance** : White Powder

**Availability** : 10 MG

**Solvent** : Ethanol

**Storage Conditions** : 4 °C

**Synonyms** : 11-Ketotestosterone  
4-Androsten-17 $\beta$ -ol-3,11-dion  
17 $\beta$ -Hydroxy-4-androstene-3,11-dion

**Other** : Standard of ELISA by 117(11-oxo-Testosterone-3-CMO, HRP) and 118E(Anti 11-Oxo-Testosterone)

## Related products

| Product Name                                                                                                                        | Maker | Cat. Num. | Quantity  |
|-------------------------------------------------------------------------------------------------------------------------------------|-------|-----------|-----------|
| Testosterone-3(E)-HRP                                                                                                               | FKA   | 101       | 1000 TEST |
| Anti Testosterone-3(E) Polyclonal Antibody                                                                                          | FKA   | 102       | 2000 TEST |
| Anti Testosterone-3(E) Polyclonal Antibody (Purified)                                                                               | FKA   | 102E      | 2000 TEST |
| Testosterone-11 $\alpha$ -HRP                                                                                                       | FKA   | 103       | 1000 TEST |
| Anti Testosterone-11 $\alpha$ Polyclonal Antibody                                                                                   | FKA   | 104       | 2000 TEST |
| Anti Testosterone-11 $\alpha$ Polyclonal Antibody (Purified)                                                                        | FKA   | 104E      | 2000 TEST |
| DHT-11 $\alpha$ -HRP                                                                                                                | FKA   | 111       | 1000 TEST |
| Anti DHT-11 $\alpha$ Polyclonal Antibody                                                                                            | FKA   | 112       | 2000 TEST |
| Anti DHT-11 $\alpha$ Polyclonal Antibody (Purified)                                                                                 | FKA   | 112E      | 2000 TEST |
| 11-Ketotestosterone (11-KT)                                                                                                         | FKA   | 117ST     | 10 MG     |
| 11-Oxo-Testosterone-3-HRP                                                                                                           | FKA   | 117       | 1000 TEST |
| Anti 11-Oxo-Testosterone-3 Polyclonal Antibody                                                                                      | FKA   | 118       | 2000 TEST |
| Anti 11-Oxo-Testosterone-3 Polyclonal Antibody (Purified)                                                                           | FKA   | 118E      | 2000 TEST |
| 11 $\beta$ -OH-Testosterone-3-HRP                                                                                                   | FKA   | 147       | 1000 TEST |
| Anti 11 $\beta$ -OH-Testosterone-3 Polyclonal Antibody (Purified)                                                                   | FKA   | 148E      | 2000 TEST |
| Anti Androstenedione-11 $\alpha$ Polyclonal Antibody                                                                                | FKA   | 106       | 2000 TEST |
| Anti Androstenedione-11 $\alpha$ Polyclonal Antibody (Purified)                                                                     | FKA   | 106E      | 2000 TEST |
| 16 $\alpha$ -OH-Androstenedione-3-HRP                                                                                               | FKA   | 123       | 1000 TEST |
| Anti 16 $\alpha$ -OH-Androstenedione-3 Polyclonal Antibody (Purified)                                                               | FKA   | 124E      | 2000 TEST |
| Androstenedione-3-HRP                                                                                                               | FKA   | 137       | 1000 TEST |
| Anti Androstenedione-3 Polyclonal Antibody                                                                                          | FKA   | 138       | 2000 TEST |
| Anti Androstenedione-3 Polyclonal Antibody (Purified)                                                                               | FKA   | 138E      | 2000 TEST |
| 11 $\beta$ -OH-Androstenedione-HRP                                                                                                  | FKA   | 139       | 1000 TEST |
| Anti 11 $\beta$ -OH-Androstenedione Polyclonal Antibody (Purified)                                                                  | FKA   | 140E      | 2000 TEST |
| 4-Androstenetrione-3-HRP                                                                                                            | FKA   | 149       | 1000 TEST |
| Anti 4-Androstenetrione-3 Polyclonal Antibody (Purified)                                                                            | FKA   | 150E      | 2000 TEST |
| Anti 5 $\alpha$ -Androstane-3 $\alpha$ , 17 $\beta$ -diol Polyclonal Antibody                                                       | FKA   | 130       | 2000 TEST |
| Anti 5 $\alpha$ -Androstane-3 $\alpha$ , 17 $\beta$ -diol Polyclonal Antibody (Purified)                                            | FKA   | 130E      | 2000 TEST |
| 5 $\alpha$ -Androstane-3 $\alpha$ , 17 $\beta$ -diol-3-Glucuronide-HRP                                                              | FKA   | 141       | 1000 TEST |
| Anti 5 $\alpha$ -Androstane-3 $\alpha$ , 17 $\beta$ -diol-3-Glucuronide Polyclonal Antibody (Purified)                              | FKA   | 142E      | 2000 TEST |
| 5 $\alpha$ -Androstane-3 $\alpha$ , 17 $\beta$ -Diol-17-Glucuronide-HRP                                                             | FKA   | 143       | 1000 TEST |
| Anti 5 $\alpha$ -Androstane-3 $\alpha$ , 17 $\beta$ -Diol-17-Glucuronide Polyclonal Antibody (Purified)                             | FKA   | 144E      | 2000 TEST |
| Anti 5 $\alpha$ -Androstane-3 $\alpha$ , 11 $\alpha$ , 17 $\beta$ -Triol-17-Glucuronide-11-Succinate Polyclonal Antibody (Purified) | FKA   | 146E      | 2000 TEST |
| DHA-11 $\alpha$ -HRP                                                                                                                | FKA   | 107       | 1000 TEST |
| Anti DHA-11 $\alpha$ Polyclonal Antibody                                                                                            | FKA   | 108       | 2000 TEST |
| Anti DHA-11 $\alpha$ Polyclonal Antibody (Purified)                                                                                 | FKA   | 108E      | 2000 TEST |
| DHA-3-HRP                                                                                                                           | FKA   | 109       | 1000 TEST |
| Anti DHA-3 Polyclonal Antibody                                                                                                      | FKA   | 110       | 2000 TEST |
| Anti DHA-3 Polyclonal Antibody (Purified)                                                                                           | FKA   | 110E      | 2000 TEST |
| 16 $\alpha$ -OH-DHA-3-HRP                                                                                                           | FKA   | 125       | 1000 TEST |
| Anti 16 $\alpha$ -OH-DHA- Polyclonal Antibody                                                                                       | FKA   | 126       | 2000 TEST |
| Anti 16 $\alpha$ -OH-DHA- Polyclonal Antibody (Purified)                                                                            | FKA   | 126E      | 2000 TEST |

*For research use only, Not for diagnostic use.*



**COSMO BIO CO., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9613

FAX : +81-3-5632-9619