

ISO Oxygen Absorbing System

Oxygen Absorber A-500HS , A-750HS , A-2500HS



This system can absorb oxygen completely and from your noble specimen such as unstable crystals, metals, medicines, important sensor probes.

You can use this system anywhere in any time promptly because this oxygen absorber is packed individually by an aluminum laminated pouch with oxygen indicator.

This absorbents absorb oxygen faster than the other one for foods use.

This system is fit for biological ,electronic ,medical use or research.

....for easily oxidizable vitamins, antibiotics, and unsaturated fatty acids.

....for natural remedies, samples of animals and plants against mold and vermin.

....for keeping quality of anaerobic medicines.

....for preventing from oxidization and rust of electrodes, medical utensils, and metal alloys.

....for transportation of specimens sensitive for atmosphere.

Table of Products

Cat.No.	Description	contents
A-500-100	Oxygen Absorber A-500-100	A-500HS 100 pcs.(with indicators)
A-500-50S	Oxygen Absorber Set	A-500HS 50 pcs.(with indicators) AP-1826 Storage Bags 25 pcs. AP-1522 Storage Bags 25 pcs. CL-21 Sealing Clips 2pcs.
AP-1826-300	Storage Bags(18×26cm)	300 pcs.
AP-1522-300	Storage Bags(15×22cm)	300 pcs.
A-750-100	Oxygen Absorber A-750-100	A-750HS 100 pcs.(with indicators)
A-2500-50	Oxygen Absorber A-2500-50	A-2500HS 50 pcs.(with indicators)
CL-21-10	Sealing Clips(21 cm)	10 pcs.

[Air volume] A-500HS: 500ml, A-750HS:750ml, A-2500HS:2500ml

Tel. 81-45-813-3851 Fax. 81-45-813-3853 E-mail: iso-co@nifty.com

I.S.O.Inc.

Rune 307, 1500 Okatsu-cho, Izumi-ku, Yokohama, 245-0003 JAPAN