

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. LSL-LB-0445

Anti type IV Collagen

Product type Primary antibodies

HostRabbitSourceSerumFormLiquid

This product does not contain preservatives such as NaN3.

Volume 100μl

Concentration

Specificity Does not cross react with Type I, II, III, V, VI collagen, laminin and fibronectin. Appears as

a single band by Western blotting.

MAP sequence match 20 out of 20 amino acids of the human and mouse.

Antigen Isotype

rat Multiple Antigen Peptide (MAP)

Application notes ELISA, WB , IHC, IF

Recommended use

According to the ELISA assay, results are positive for dilutions up to 500,000 fold against

linear synthetic peptides.

Immunohistochemistry on Frozen or Paraffin section

Recommended dilutions

Western Blotting; more than 1/1,000

Immunofluoresence; more than 1/500 (recommended: 1/200~1/2,000)

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity The MAP contains the SCL sequence of type IV Collagen (α1 IV) MAP sequence match

18/20 aa \sim 13/20 aa of the human, mouse and bovine.

Cross reaction will differ between species.

Storage Store below -20°C (below -70°C for prolonged storage).

After thawing, store in small aliquots in sealable vials and store below -70°C. To prevent degradation from repeated thawing, store the antiserum between 0 to 4°C after second

thawing.

Stable for three years at -70°C.

References

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp e-mail: export@cosmobio.co.jp

[Outside Japan] Phone: +81-3-5632-9617 [国内連絡先] Phone: +81-3-5632-9610 FAX: +81-3-5632-9618 FAX: +81-3-5632-9619

Anti (MAP) type N Collagen

I. 内容

Lot No. 738030

免	疫	動	物	ウサギ
性状・包装サイス				全血清・100μℓ
力			価	ELISA で500,000 倍希釈まで陽性(合成直鎖ペプチ ドに対して)
抗	原	由	来	Multiple Antigen Peptide *
種	間	交	差	ヒト、サル、ラット、マウス、ウシ等と交差*
特			徴	Laminin,type I, II, III, V, VI collagen & Fibronectinとは非交差。 Western Blotting で type IV Collagen (human & mouse)の分子量に相当する位置に陽性バンドを確認。免疫組織染色で基底膜と良く反応。Protease処理でその反応性は著しく向上する。又、Native Collagen に対する反応性は低い。
標準	基 希系	平		蛍光抗体法で1:500以上(1/200 ~1/2,000) Western Blottingに関しては1:1,000以上

()内は推薦希釈率

*種によって交差性が異なると思われます、希釈率等充分に検討下さい。但し、MAPに合成した type IV Collagen (α 1(IV))中の-SCL-を含むアミノ酸配列(20残基)はヒト、マウスで20/20一致している。又、 α 2 \sim 6(IV) に関してはヒト、マウス、ウシなどで18/20 \sim 13/20一致している。

Ⅱ. 保存上の注意

-20℃以下(長期間の場合は-70℃以下)で凍結して下さい。 解凍後は密栓のできる小型容器に研究の規模に応じて少量づつ分注し、-70℃以下で保存して下さい。 凍結融解の繰り返しによる力価の低下を避ける為、再解凍後の抗血清 は0~4℃に保ち操作・保存して下さい。

Ⅲ. 安定性

-70℃で3年間安定。 但し、NaN。等の防腐剤は入っていません。

製造元

総発売元

\$1 株式会社 エル・エス・エル

コスモ・バイオ株式会社