

Rep-HepG2 Maintenance Medium

【Code No. REP-HEPG2M】

Feb 7, 2025

※For research use only. Not for clinical diagnosis.

《 I . Product Overview 》

Rep-HepG2 Maintenance Medium optimized for culturing Rep-HepG2 cells 【Catalog No.: REP-HEPG2C】 .

Components

Components	size	Quantity	Storage
Rep-HepG2 Maintenance Medium	125 mL	1	－ 20℃
supplements	125 µL	1	－ 20℃

《 II . Storage Conditions 》

All frozen products are shipped in dry ice. When you receive the shipment, make sure to store at -20℃immediately.

Thaw the medium at 4 °C before use

Supplement should be thawed just before use.

《 III . Medium preparation 》

Add 125ul of supplement to 125ml of the media.

Prepared media can be stored at 4℃ for up to 1 to 2 weeks.

If prepared media is not to be used immediately, it is acceptable to store aliquots of it at -20℃.

※Avoid repeated freeze/thaw cycles.

Rep-HepG2 Maintenance Medium

【Code No. REP-HEPG2M】

2025 年 2 月 7 日作成

※本マニュアルをご精読のうえ、研究目的にのみご使用ください。

《 I . 製品情報 》

本製品は Rep-HepG2 細胞 (品番 : REP-HEPG2C) 培養用の専用メディウムです。

製品構成

構成	容量	本数	保存温度
Rep-HepG2 Maintenance Medium	125 mL	1 本	−20℃
サプリメント	125 μ L	1 本	−20℃

《 II . 保管方法 》

本製品は、冷凍状態で発送しております。本製品受領後は速やかに−20℃以下で保管してください。
培地は使用前にあらかじめ4℃で解凍してください。
サプリメントは使用直前に解凍してください。

《 III . 培地調製方法 》

メディウム1本(125mL)に対しサプリメント1本(125 μ L)を添加してください。
サプリメント添加後の培地は4℃で1〜2週間保管可能です。
一度に使用しない場合には、サプリメント添加後ただちに小分けして凍結保存してください。
※培地温めの繰り返し、凍結融解の繰り返しは避けてください。