## rTaq (recombinant Taq) DNA Polymerase



Code No. TAP-211

Lot No. \*\*\*\*\*

Storage Store at -20°C

Size 250units

Source : Escherichia coli JM109(pTQ1037)

Concentration :  $\underline{\phantom{a}}$  units/ $\mu$ l

Unit Definition : One unit of enzyme is defined as the amount of enzyme that will incorporate

10 nmoles of dNTPs into acid insoluble material in 30 minutes at 75°C.

Assay Condition : 25 mM TAPS(pH9.3)

50 mM KCl 2 mM MgCl<sub>2</sub>

200 μM each of dATP,dGTP,dTTP

100 μM [ $\alpha$ -<sup>32</sup>P]-dCTP

20 μg activated salmon sperm DNA per 50 μl reaction

Storage Buffer : 20 mM Tris-HCl(pH8.0)

100 mM KCI 0.1 **EDTA** mM DTT 1 mM Tween-20 0.5 0.5 % Nonidet P-40 50 Glycerol

rTaq DNA Polymerase : 500 mM KCl

 $10 \times Buffer(+Mg)$  100 mM Tris-HCl(pH8.3 at 25°C)

15 mM MgCl<sub>2</sub>

dNTPs : 2 mM dATP,dGTP,dCTP,dTTP each

**Quality Control** 

1. Endonuclease Activity : When 40 units of this enzyme were incubated with 1  $\mu$ g of  $\lambda$  -DNA for 16 hours

at 75°C, no endonuclease activity was observed after agarose gel electrophoresis.

2. Nicking Activity : When 40 units of this enzyme were incubated with 1 µg of pBR322 for 16 hours

at 75°C, no nicking activity was observed after agarose gel electrophoresis.

Purchase of this product includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research. No other patent rights (such as 5' Nuclease Process patent rights) are conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.



## COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618 / e-mail: export@cosmobio.co.jp / URL:www.cosmobio.com