



## C-Slide®

Count four samples at once - no washing and no risk of contamination  
Designed for precise results and greater efficiency

## ■ C-Slide® Key Feature

### Efficient

No need to wash the hemocytometer after cell counting, which saves time and effort.

### Accurate

The consistent depth of the C-slide ensures exact volume control, leading to accurate cell counting.

### Affordable

The 4-channel design enables a larger number of experiments and eliminates the need for coverslips.

### Safe

Reduces the risk of exposure to infectious or harmful samples.



Individual package | 50 slides / box (200 tests)



## ■ Specification

Dimension	25 × 75 × 2 mm
Grid	Neubauer Improved (NI)
Loading volume	20 µL / chamber
Packing unit	50 slides / box
Total tests	200 tests / box

Curiosis Inc.

+82 2 508 5237 | sales@curiosis.com | South Korea

FOR RESEARCH USE ONLY and not for use in diagnostic procedures. Specifications subject to change without notice.



CRB001 -2304Rev01