

Di-ubiquitin (K27-linked) [untagged]

Ubiquitin/Ubiquitin-Like Protein Dimer



Cat. No. 60-0103-010

Lot. No. 30084

Quantity: 10 µg

Storage: -70°C

FOR RESEARCH USE ONLY

NOT FOR USE IN HUMANS

CERTIFICATE OF ANALYSIS Page 1 of 2

Background

Ubiquitin (Ub) is a highly conserved 76 amino-acid protein found throughout eukaryotic cells. A vast number of cellular processes, including targeted protein degradation, cell cycle progression, DNA repair, protein trafficking, inflammatory response, virus budding, and receptor endocytosis, are regulated by Ub-mediated signalling; where the target protein is tagged by single or multi-monomeric Ub (monomeric Ub attached to multiple sites on the substrate) or a polymeric chain of Ubs (Fushman *et al.*, 2010). This post-translational modification is tightly controlled by an enzymatic cascade involving several enzymes (E1, E2, and E3) and occurs through either an isopeptide bond between the C-terminal Glycyl residue of Ub and the epsilon amino group of a Lysyl residue on a target protein or through a peptide bond between the C-terminal Glycyl residue of Ub and the N-terminal amine on a further Ub. In the former (isopeptide bond-linked) case the substrate protein may either be ubiquitin itself – thus leading to the generation of poly-ubiquitin chains – or another target protein (Fushman *et al.*, 2010). Thus, ubiquitin can be attached to a substrate either as a monomer or as a poly-ubiquitin chain. Further – depending on their linkage type (M1, K6, K11, K27, K29, K33, K48 and K63 linked) – the Ub chains can take different structural forms. Chains containing all eight possible Ub linkages have been found in living cells and different ubiquitin chain types may encode different biological signals, allowing this single protein to mediate many diverse functions (Komander 2009; Weeks *et al.*, 2009; Walczak *et al.*, 2012). The functionality of Ub chains is most commonly associated with their attachment to substrate proteins but there is also evidence that they may also play a role in cellular signalling as free chains (Braten *et al.*, 2012).

Continued on page 2

Physical Characteristics

Protein Sequence:

```
MQIFVKTLTGKTITLEVEPSDTIENVKAKIQDKEGIPPDQORLIFAGKQLEDGRTLSDYNIQKESTLHLVLRGG  
MQIFVKTLTGKTITLEVEPSDTIENVKAKIQDKEGIPPDQORLIFAGKQLEDGRTLSDYNIQKESTLHLVLRGG  
K27
```

Species: human

Molecular Weight: 17.1 kDa

Source: synthetic/chemical ligation

Purity: >98% by InstantBlue™ SDS-PAGE

Quantity: 10 µg

Stability/Storage: 12 months at -70°C; aliquot as required

Concentration: 0.5 mg/ml


Accession Number: P62987

Formulation: 50 mM HEPES pH 7.5, 150 mM NaCl₂, 2 mM DTT, 10% Glycerol

Quality Assurance

Purity:

4-12% gradient SDS-PAGE
InstantBlue™ staining
Lane 1: MW markers
Lane 2: 1 µg Di-ubiquitin (K27-linked)




Purity of the linkage type:

The linkage type (K27) was confirmed by tandem mass spectrometry.

Di-ubiquitin cleavage assay:

The capacity of the di-ubiquitin substrate to be cleaved was tested using a promiscuous – with respect to ubiquitin linkage specificity – deubiquitylase (GST-USP2). Incubation of the di-ubiquitin for 1 hour at 37°C was compared either in the absence (Lane 2) or presence (Lane 3) of GST-USP2. The reaction products were compared alongside two control samples containing either mono-ubiquitin (Lane 4) or GST-USP2 (Lane 5) only. Cleavage of the di-ubiquitin and generation of mono-ubiquitin was determined by running reactions on a 4-12% SDS-PAGE gel and staining with InstantBlue™ (Lane 1; molecular weight markers).



www.ubiquigent.com
Dundee, Scotland, UK

ORDERS / SALES SUPPORT

International: +1-617-245-0003
US Toll-Free: 1-888-4E1E2E3 (1-888-431-3233)
Email: sales.support@ubiquigent.com

UK HQ and TECHNICAL SUPPORT

International: +44 (0) 1382 381147 (9AM-5PM UTC)
US/Canada: +1-617-245-0020 (9AM-5PM UTC)
Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2012. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0

Di-ubiquitin (K27-linked) [untagged]

Ubiquitin/Ubiquitin-Like Protein Dimer



Cat. No. 60-0103-010

Lot. No. 30084

Quantity: 10 µg

Storage: -70°C

FOR RESEARCH USE ONLY

NOT FOR USE IN HUMANS

CERTIFICATE OF ANALYSIS Page 2 of 2

Background

Continued from page 1

A mass spectrometry-based study found that K27 linkages account for 9% of all yeast ubiquitin-ubiquitin linkages. The relative abundance of the other linkages were K6 (11%), K11 (28%), K29 (3%), K33 (4%), K48 (29%) and K63 (16%) (Xu *et al.*, 2009). The role of the K27 linkage remains elusive, possibly in part due to the close proximity of K29 and K33 linkages making MS analysis a challenge (owing to the nature of the tryptic digest fragments). No cellular role has been associated with K27-linked polymers to date (Komander 2009).

References:

Braten O, Shabek N, Kravtsova-Ivantsiv Y, Ciechanover A (2012) Generation of free ubiquitin chains is upregulated in stress, and facilitated by the HECT domain ubiquitin ligases UFD4 and HUL5. *Biochem J* **444**, 611-617.

Fushman D, Walker O (2010) Exploring the linkage dependence of polyubiquitin conformations using molecular modeling. *Journal of Molecular Biology* **395**, 803-814.

Komander D (2009) The emerging complexity of protein ubiquitination. *Biochem Soc Trans* **37**, 937-953.

Walczak H, Iwai K, Dikic I (2012) Generation and physiological roles of linear ubiquitin chains. *BMC Biol* **10**, 23.

Weeks SD, Grasty KC, Hernandez-Cuebas L, Loll PJ (2009) Crystal structures of Lys-63-linked tri- and di-ubiquitin reveal a highly extended chain architecture. *Proteins* **77**, 753-759.

Xu P, Duong DM, Seyfried NT, Cheng D, Xie Y, *et al.*, (2009) Quantitative proteomics reveals the function of unconventional ubiquitin chains in proteasomal degradation. *Cell* **137**, 133-145.



www.ubiquigent.com
Dundee, Scotland, UK

ORDERS / SALES SUPPORT

International: +1-617-245-0003
US Toll-Free: 1-888-4E1E2E3 (1-888-431-3233)
Email: sales.support@ubiquigent.com

UK HQ and TECHNICAL SUPPORT

International: +44 (0) 1382 381147 (9AM-5PM UTC)
US/Canada: +1-617-245-0020 (9AM-5PM UTC)
Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2012. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0