# **UBE2T** (HSPC150) [untagged]

E2 – Ubiquitin Conjugating Enzyme

Alternate Name: HSPC150

Cat. No. 62-0070-020 Lot. No. 30195

FOR RESEARCH USE ONLY

Quantity: 20 µg -70°C Storage:

NOT FOR USE IN HUMANS



**CERTIFICATE OF ANALYSIS Page 1 of 2** 

## **Background**

The enzymes of the ubiquitylation pathway play a pivotal role in a number of cellular processes including regulated and targeted proteasomal degradation of substrate proteins. Three classes of enzymes are involved in the process of ubiquitylation; activating enzymes (E1s), conjugating enzymes (E2s) and protein ligases (E3s). UBE2T is a member of the E2 conjugating enzyme family and cloning of the human gene was first described by Zhang et al. (2000). UBE2T is integral to the Fanconi Anemia pathway for DNA damage repair. UBE2T binds to the C-terminal PH domain of FAN-CL the ubiquitin ligase subunit of the Fanconi Anemia (FA) core complex, which leads to the monoubiquitylation of FANCD2 and FANCI (Longerich et al., 2009; Machida et al., 2006). E3 ligase activity is not determined by assembly of the FA core complex but by the DNA damage-induced subcellular localization of the complex to chromatin. UBE2T and FANCD2 access this subcellular fraction independently and FANCD2 monoubiquitylation is regulated by the formation of an E2/E3 holoenzyme on chromatin. DNA damage in UBE2T-depleted human osteosarcoma cells leads to the formation of abnormal chromosomes that are a hallmark of FA (Alpi et al., 2007). UBE2T expression has been analysed in lung cancer tissue and compared to normal human tissue.

Continued on page 2

## **Physical Characteristics**

Species: human

Source: E. coli expression

Quantity: 20 µg

Concentration: 1 mg/ml

Formulation: 50 mM HEPES pH 7.5. 150 mM sodium chloride, 2 mM dithiothreitol, 10% glycerol

Molecular Weight: ~23 kDa

Purity: >80% by InstantBlue™ SDS-PAGE

Stability/Storage: 12 months at -70°C;

aliquot as required

#### **Protein Sequence:**

 ${\tt GPGS} \textbf{\textit{M}} {\tt QRASRLKRELHMLATEPPPGITCWQD}$ KDQMDDLRAQILGGANTPYEKGVFKLEVI IPERYPFEPPQIRFLTPIYHPNIDSAGRI CLDVLKLPPKGAWRPSLNIATVLTSIQLLM SEPNPDDPLMADISSEFKYNKPAFLKNARQW TEKHARQKQKADEEEMLDNLPEAGDSRVHN STOKRKASQLVGIEKKFHPDV

The residues underlined remain after cleavage and removal of the purification tag


UBE2T (regular text): Start bold italics (amino acid residues

Accession number: NP\_054895

# **Quality Assurance**

#### **Purity:**

4-12% gradient SDS-PAGE InstantBlue™ staining Lane 1: MW markers Lane 2: 1 µg UBE2T



### **Protein Identification:**

Confirmed by mass spectrometry.

#### **E2-Ubiquitin Thioester Loading Assay:**

The activity of UBE2T was validated by loading E1 UBE1 activated ubiquitin onto the active cysteine of the UBE2T E2 enzyme via a transthiolation reaction. Incubation of the UBE1 and UBE2T enzymes in the presence of ubiquitin and ATP at 30°C was compared at two time points, T<sub>0</sub> and T<sub>10</sub> minutes. Sensitivity of the ubiquitin/UBE2T thioester bond to the reducing agent DTT was confirmed.



Dundee, Scotland, UK

#### **ORDERS / SALES SUPPORT**

US Toll-Free: 1-888-4E1E2E3 (1-888-431-3233) International: +1-617-245-0003

Email: sales.support@ubiquigent.com

#### **UK HQ and TECHNICAL SUPPORT**

International: +1-617-245-0020 Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© **Ubiquigent 2013**. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0

# **UBE2T** (HSPC150) [untagged]

E2 - Ubiquitin Conjugating Enzyme

Alternate Name: HSPC150

**Cat. No. 62-0070-020** Quantity: **Lot. No. 30195** Storage:

FOR RESEARCH USE ONLY NOT FOR USE IN HUMANS

20 µg

-70°C



**CERTIFICATE OF ANALYSIS Page 2 of 2** 

### **Background**

Continued from page 1

UBE2T was found to be significantly upregulated at both the protein and mRNA level suggesting involvement in the malignant cell phenotype (Hao *et al.*, 2008).

#### References:

Alpi A, Langevin F, Mosedale G, Machida YJ, Dutta A, Patel KJ (2007) UBE2T, the Fanconi anemia core complex, and FANCD2 are recruited independently to chromatin: a basis for the regulation of FANCD2 monoubiquitination. *Mol Cell Biol* 27, 8421-30.

Hao J, Xu A, Xie X, Tian T, Gao S, Xiao X, He D (2008) Elevated expression of UBE2T in lung cancer tumors and cell lines. *Tumour Biol* **29**, 195-203.

Longerich S, San Filippo J, Liu D, Sung P (2009) FANCI binds branched DNA and is monoubiquitinated by UBE2T-FANCL. *J Biol Chem* **284**, 23182-6.

Machida YJ, Machida Y, Chen Y, Gurtan AM, Kupfer GM, D'Andrea AD, Dutta A (2006) UBE2T is the E2 in the Fanconi anemia pathway and undergoes negative autoregulation. *Mol Cell* 23, 589-96

Zhang QH, Ye M, Wu XY, Ren SX, Zhao M, Zhao CJ, Fu G, Shen Y, Fan HY, Lu G, Zhong M, Xu XR, Han ZG, Zhang JW, Tao J, Huang QH, Zhou J, Hu GX, Gu J, Chen SJ, Chen Z (2000) Cloning and functional analysis of cDNAs with open reading frames for 300 previously undefined genes expressed in CD34+ hematopoietic stem/progenitor cells. *Genome Res* 10, 1546-60.



Dundee, Scotland, UK

#### **ORDERS / SALES SUPPORT**

**US Toll-Free:** 1-888-4E1E2E3 (1-888-431-3233) **International:** +1-617-245-0003

Email: sales.support@ubiquigent.com

#### **UK HQ and TECHNICAL SUPPORT**

International: +1-617-245-0020

Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2013. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0