

USP7 [6His-tagged]

Deubiquitylating Enzyme

Alternate Names: Herpesvirus-Associated Ubiquitin-Specific Protease, HAUSP VMW110-associated protein

Cat. No. 64-0003-050
Lot. No. 1736

Quantity: 50 µg
Storage: -70°C

FOR RESEARCH USE ONLY

NOT FOR USE IN HUMANS



CERTIFICATE OF ANALYSIS Page 1 of 2

Background

The deubiquitylating enzymes (DUBs) regulate ubiquitin dependent signaling pathways. The activities of the DUBs include the generation of free ubiquitin from precursor molecules, the recycling of ubiquitin following substrate degradation to maintain cellular ubiquitin homeostasis and the removal of ubiquitin or ubiquitin-like proteins (UBL) modifications through chain editing to rescue proteins from proteasomal degradation or to influence cell signalling events (Komander *et al.*, 2009). There are two main classes of DUB, cysteine proteases and metalloproteases. Ubiquitin specific processing protease 7 (USP-7) is a member of the cysteine protease enzyme family and cloning of the gene in humans was first described by Everett *et al.* (1997). Overexpression of p53 and USP7 stabilizes p53 through the removal of ubiquitin moieties from polyubiquitylated p53 (Kon *et al.*, 2010). Inhibition of USP7 expression results in the accumulation and stabilisation of p53 protein. However, USP7 has been shown to stabilise mdm2 a ubiquitin ligase that promotes the degradation of p53 (Holowaty and Frappier. 2004; Krishna and Grishin. 2004). Thus USP7 appears to play multiple roles in regulating the p53/mdm2 pathway and maintaining steady-state levels of p53 in the cell. Due to the role USP7 has been shown to play in the regulation of p53 its importance as a thera-

Continued on page 2

Physical Characteristics

Species: human

Source: *E. coli* expression

Quantity: 50 µg

Concentration: 0.5 mg/ml

Formulation: 50 mM HEPES pH 7.5,
150 mM sodium chloride, 2 mM
dithiothreitol, 10% glycerol

Molecular Weight: ~130 kDa


Purity: >80% by InstantBlue™ SDS-PAGE

Stability/Storage: 12 months at -70°C;
aliquot as required

Protein Sequence: Please see page 2

Quality Assurance

Purity:
4-12% gradient SDS-PAGE
InstantBlue™ staining
Lane 1: MW markers
Lane 2: 1 µg His-USP7



Protein Identification:

Confirmed by mass spectrometry.

Deubiquitylating Enzyme Assay:

The activity of His-USP7 was validated by determining the increase in fluorescence measured as a result of the enzyme catalysed cleavage of the fluorogenic substrate Ubiquitin-Rhodamine110-Glycine generating Ubiquitin and Rhodamine110-Glycine. Incubation of the substrate in the presence or absence of His-USP7 was compared confirming the deubiquitylating activity of His-USP7.



www.ubiquigent.com
Dundee, Scotland, UK

ORDERS / SALES SUPPORT

International: +1-617-245-0003
US Toll-Free: 1-888-4E1E2E3 (1-888-431-3233)
Email: sales.support@ubiquigent.com

UK HQ and TECHNICAL SUPPORT

International: +44 (0) 1382 381147 (9AM-5PM UTC)
US/Canada: +1-617-245-0020 (9AM-5PM UTC)
Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2011. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0

USP7 [6His-tagged]

Deubiquitylating Enzyme

Alternate Names: Herpesvirus-Associated Ubiquitin-Specific Protease, HAUSP VMW110-associated protein

Cat. No. 64-0003-050

Lot. No. 1736

Quantity: 50 µg

Storage: -70°C



FOR RESEARCH USE ONLY

NOT FOR USE IN HUMANS

CERTIFICATE OF ANALYSIS Page 2 of 2

Background

Continued from page 1

peutic target for treating hematopoietic tumours has been highlighted (Cheon and Baek, 2006).

References:

Cheon KW, Baek KH (2006) HAUSP as a therapeutic target for hematopoietic tumors (review). *Int J Oncol* **28**, 1209-15.

Everett RD, Meredith M, Orr A, Cross A, Kathoria M, Parkinson J (1997) A novel ubiquitin-specific protease is dynamically associated with the PML nuclear domain and binds to a herpesvirus regulatory protein. *EMBO J* **16**, 1519-30.

Holowaty MN, Frappier L (2004) HAUSP/USP7 as an Epstein-Barr virus target. *Biochem Soc Trans* **32**, 731-2.

Komander D, Clague MJ, Urbe S (2009) Breaking the chains: structure and function of the deubiquitinases. *Nat Rev Mol Cell Biol* **10**, 550-63.

Kon N, Kobayashi Y, Li M, Brooks CL, Ludwig T, Gu W (2010) Inactivation of HAUSP in vivo modulates p53 function. *Oncogene* **29**, 1270-9.

Krishna SS, Grishin NV (2004) The finger domain of the human deubiquitinating enzyme HAUSP is a zinc ribbon. *Cell Cycle* **3**, 1046-9.

Physical Characteristics

Continued from page 1

Protein Sequence:

M G S S H H H H H S S G L E V L F Q G P G S M N
HQQQQQQQKAGEQQQLSEPEDMEMEAGDTD
DPPRITQNPVINGNVALSDGHNTAEEDMED
DTSWRSEATFQFTVERFSRLSESVLSPP
CFVRNLPWKIMVMPrfYPDRPHQKSVGF
FLQCNAESDSTSWSCHAQAVLKIINyRD
DEKSFsRRISHLFFHKENDWGFsNFMawSEVT
DPEKGFIDDDKVTfEVfVQADAPHGVAWDSK
KHTGYVGLKNGGATCYMNSLLQLFFTNQL
RKAVYMPTEGDDSSKSVPLALQRFVYELQHS
DKPVGTKKLTkSfGWETLDSFMQHDVQEL
CRVLLDNVENKMGTCVEGTIPKLFrgKM
VSYIQCKEVDYrSDRREDYDIQLSikGK
KNIFESFVDYVAVEQLDGNKYDAGEHGLQE
AEKGVKFLTLPPVLHLQLMRfMYDPQTDQNI
KINDRFEfPEQLPLDEFLQKTDPKDPANYIL
HAVLVHSGDNHGGHYVVYLNPkGDGKW
CKFDDDVVSrCTKEEAIEHNYGGHDDDLs
VRHCTNAYMLVYIRESKLSEVLQAVTDHDI
PQQLVERLQEEKRIEaQKRKERQEAHLYMQVQ
IVAEDQFCGHQGNdMYDEEKVKYtVfKVLKNS
SLAEfVQSLsQTMGFpQDQIRLWPMQARSNGT
KRPAmlDNEADGNkTMIeLSDNENPWTIFLET
VDPELAASGATLPKfDKDHDVMLfLkMYDPK
TRSLNYCGHIYtPIsCKIRDLLPVMCDRAG
FIQDTSLILYEEVKPNLTERIQDYDVSLDKA
LDELMDGDIIVfQKDDPENDNSelPTakeY
FRDLYHRVDVIFCDKTIpNDPGfVVTLsNRM
NYfQVAKtVAQRlNTDpMLLQfFKSQGYRDG
PGNPLRHNYEGTLRDLLQfFKPRQPKKLYY
QQLKMKITDFENRRSfKCIWLNsqFREEEIT
LYPDKHGcVRDLLEECkKAVELGekASGKL
RLLEIVSYKIIGVHQEDELLECLSPATSRT
FRIEEIPLDQVDIDKENEMLVtVAHFHKEV
FGTfGIPfLLRIHQGEHFREVMKRIQSLD
DIQEKEfEKfKFAIVMTGRHQYINeDEYEVN
LKDFEPQPGNMShPRPWLGLDHFNKAPKRSRY
TYLEKAIKIHN*

Tag (bold text): N-terminal His

Protease cleavage site: PreScission™ (LEVLfQ▼GP)

USP7 (regular text): Start **bold italics** (amino acid residues 1-1102)

Accession number: CAA96580



www.ubiquigent.com
Dundee, Scotland, UK

ORDERS / SALES SUPPORT

International: +1-617-245-0003
US Toll-Free: 1-888-4E1E2E3 (1-888-431-3233)
Email: sales.support@ubiquigent.com

UK HQ and TECHNICAL SUPPORT

International: +44 (0) 1382 381147 (9AM-5PM UTC)
US/Canada: +1-617-245-0020 (9AM-5PM UTC)
Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2011. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0