



**RANDOX** *life sciences*

## PRODUCT SPECIFICATION

<b>Product</b>	:	Anti-Creatinine
<b>Cat. No.</b>	:	PAS 9686
<b>Lot No.</b>	:	RS300-24
<b>Immunogen</b>	:	Creatinine(N)-BSA
<b>Host Animal</b>	:	Sheep
<b>Format</b>	:	Ig fraction
<b>Ig Concentration</b>	:	6.13mg/ml (U.V. abs @ 280nm)
<b>Working Concentration</b>	:	2µg/ml (EIA)
<b>PH</b>	:	7.2
<b>Buffer</b>	:	20mM Phosphate, 150mM Sodium Chloride
<b>Preservative</b>	:	0.09% Sodium Azide
<b>Storage</b>	:	-20°C

### **Randox Life Sciences**

55 Diamond Road, Crumlin, Co. Antrim, United Kingdom, BT29 4QY  
t: +44 (0) 289442 2413 f: +44 (0) 289445 2912 e: lifesciences@randox.com  
[www.randox-lifesciences.com](http://www.randox-lifesciences.com)