



Cyclin-dependent kinase 4

Human, Recombinant, *E. coli*

Cat. No.: RCP9234
RCP9235
RCP9111

Size: 10µg
100µg
1000µg

Synonym: CDK4


Description: Cyclin-dependent kinase 4 combines with the protein cyclin D to form a protein complex that promotes the passage of cells through the G1 checkpoint in the cell growth cycle. This cellular activity is regulated by the protein p16 produced by p16^{INK4a}. Normally, p16 controls cell growth by inhibiting the activity of the CDK4–cyclin D complex and stopping cells at the G1 checkpoint.

RANDOX recombinant human CDK4 is expressed in *E. coli* with an amino-terminal hexahistidine tag. This product is for research use only and is not intended for diagnostic or therapeutic use.

Form: Liquid.

Molecular Weight: 38 kDa

Purity: >95% by SDS-PAGE



References: Wölfel T. *et al.* (1995) *Science* **269**: 1281-84.