



Inhibin beta A Subunit Human, Recombinant, *E. coli*

Cat. No.: RCP9030 **Size:** 10µg

Synonym: Inhibin beta-1, Follicle-Stimulating Hormone Releasing Protein; Follicle-Stimulating Hormone Releasing Protein; Follicle-Stimulating Hormone Releasing Factor; Erythroid Differentiation Factor; EDF, Activin beta-A chain.

Description: Inhibins are dimeric peptide hormones produced by female ovarian granulose cells and male Sertoli cells as well as a variety of other tissues. Inhibins have two isoforms, A and B, with the same alpha subunit but different beta subunits. Inhibin A is a dimer of alpha and beta A subunits, inhibin B is a dimer of alpha and beta B subunits. Inhibins are thought to inhibit the production of follicle-stimulating hormone (FSH) by the pituitary gland. In addition, inhibins are also thought to play a role in the control of gametogenesis, and embryonic and foetal development.

RANDOX recombinant Inhibin beta A subunit comprises a 116 amino acid fragment (311-426) corresponding to the mature Inhibin beta A subunit protein and is expressed in *E. coli* with an amino-terminal hexahistidine tag. This product is for research use only and is not intended for diagnostic or therapeutic use.

Form: Liquid.
Supplied in 1x Laemmli Buffer (25mM Tris-HCl pH6.8, 50mM DTT, 1% (w/v) SDS, 0.1% (w/v) Bromophenol Blue, 2.5% Glycerol).

Purity: Single band on Western blot.

References: Burger, H.G. (1988) *J. Endocrinol.* **117(2)** : 159-160.
Tanimoto, K. et al., (1991) *DNA Seq.* **2** : 103-110.

Randox Life Sciences

55 Diamond Road, Crumlin, Co. Antrim, United Kingdom, BT29 4QY
t: +44 (0) 289442 2413 f: +44 (0) 289445 2912 e: lifesciences@randox.com
www.randox-lifesciences.com