



PAPP-A Fragment

Human, Recombinant, *E. coli*

Cat. No.:	RCP9037	Size:	10µg
	RCP9038		100µg
	RCP9145		1000µg


Synonym: Pappalysin-1; Pregnancy-Associated Plasma Protein-1; Insulin-like growth factor-dependent IGF binding protein-4, EC 3.4.24.79; IGF-dependent IGFBP-4 protease; IGFBP-4ase.

Description: PAPP-A is a large zinc binding protein, which acts as a metalloprotease and has been shown to cleave insulin-like growth factor binding protein-4. PAPP-A can act as a regulator of IGF bioactivity in several biological systems, including the human ovary and cardiovascular systems. Recent studies have shown that levels are elevated in patients with unstable angina or acute myocardial infarction. PAPP-A is also produced in high concentrations during pregnancy and is released into the maternal circulation. Low levels of PAPP-A have been linked with a number of foetal chromosomal abnormalities, as well as pre-eclampsia and stillbirth.

RANDOX recombinant PAPP-A comprises a 191 amino acid fragment (81-271) corresponding to the PAPP-A 'Jelly-Roll' domain fragment and is expressed in *E. coli* with an amino-terminal hexahistidine tag. This product is for research use only and is not intended for diagnostic or therapeutic use.

Form: Liquid

Purity: >95% by SDS-PAGE



References: Conover, C.A. *et al.*, 2004 *Development* **131** (5) 1187-1194