



**S100A10**

**Human, Recombinant, *E. coli***

**Cat. No.:** RCP9219  
RCP9220  
RCP9106

**Size:** 10µg  
100µg  
1000µg

**Synonym:** Calpactin light chain, p11, CAL 12, CLP11, p10, Ca [1], 42C


**Description:** The S100 family of proteins contain two EF-hand calcium-binding motifs and are thought to be involved in the regulation of a number of cellular processes including cell cycle progression and differentiation.

**RANDOX** recombinant human S100A10 is expressed in *E. coli* with an amino-terminal hexa-histidine tag. This product is for research use only and is not intended for diagnostic or therapeutic use.

**Form:** Liquid.

**Molecular Weight:** 15.3kDa

**Purity:** >95% by SDS-PAGE



**References:** Kondo A. *et al.* (2002) *Acta Med.* **56(1)**: 31-34.