



**COSMO BIO CO., LTD.**  
Inspiration for Life Science

# Antibodies

Validated for Immunohistochemistry

**Updated August 2009**

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
KAL-KR099	Anti 17βHSD type8	25 UG	Human	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KH043	Anti 3-DG-Imidazolone	50 UG	-	Mouse	JNH-27	IHC	monoclonal	IgG1	Affinity Purified
KAL-KE020	Anti 4F2 Heavy Chain (4F2hc : CD98)	250 UG	Human Rat	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
KAL-KE028	Anti 4F2 Heavy Chain (4F2hc : CD98)	25 UG	Rat	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
NOF-N213730-EX	Anti 4-hydroxy-2-hexenal	30 UG	-	Mouse	HHE53	IHC	monoclonal	IgG1k	Ig fraction - Ammonium Sulphate
CAC-NM-DND-002	Anti 6-4PPs	1 VIAL		Mouse	64M-2	ELISA IC IHC	monoclonal	IgG2ak	Semi Purified
NOF-N213810-EX	Anti 7-ketocholesterol	100 UG	-	Mouse		IHC(f)	monoclonal	IgG	Ig fraction - Protein A
NOF-N213820-EX	Anti 7-ketocholesterol	20 UG	-	Mouse		IHC(f)	monoclonal	IgG	Ig fraction - Protein A
YII-Y260-EX	Anti Acidic FGF (1-15)	50 UL	Bovine	Rabbit		IHC RIA	polyclonal		Serum
NOF-N213310-EX	Anti Acrolein	100 UG	-	Mouse		IHC	monoclonal	IgG1k	Ig fraction - Protein A
NOF-N213320-EX	Anti Acrolein	20 UG	-	Mouse		IHC	monoclonal	IgG1k	Ig fraction - Protein A
YII-Y350-EX	Anti ACTH	50 UL	Rat	Rabbit		IHC	polyclonal		Serum
YII-Y352-EX	Anti ACTH (1-23)	50 UL	Mouse Rat	Rabbit		IHC(p) IHC(f) RIA	polyclonal		Serum
YII-Y351-EX	Anti ACTH (21-39)	50 UL	Mouse Rat	Rabbit		IHC(p) IHC(f) RIA	polyclonal		Serum
KAL-KR057	Anti Addicisin	25 UG	Mouse	Rabbit		WB IHC(f)	polyclonal	IgG	Affinity Purified
KAL-KW252	Anti AFP	100 UG	Human	Mouse		IHC(p)	monoclonal	IgG2a	Ig fraction - Protein G
KAL-KH001	Anti AGEs	10 UG	-	Mouse	6D12	IHC ELISA	monoclonal	IgG1	Affinity Purified
KAL-KH001-01	Anti AGEs Biotin	10 UG	-	Mouse	6D12	IHC ELISA	monoclonal	IgG1	Affinity Purified
KAL-KH001-04	Anti AGEs FITC	10 UG	-	Mouse	6D12	IHC ELISA	monoclonal	IgG1	Affinity Purified
KAL-KH001-02	Anti AGEs HRP	20 UG		Mouse	6D12	WB IHC ELISA	monoclonal	IgG1	Ig fraction - Protein G

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
HKD-AB01	Anti Amelogenin	100 UL	Mouse Rat Bovine	Rabbit		WB IHC	polyclonal	-	Serum
YII-YA010-EX	Anti Amyloid $\beta$ Protein (1-40)	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
KAL-KW149	Anti Annexin IV	100 UG	Human Mouse Rat	Rabbit		WB IHC(p) IC	polyclonal	IgG	Affinity Purified
KAL-KW150	Anti Annexin V	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p) IC	polyclonal	IgG	Affinity Purified
YII-Y330-EX	Anti ANP	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
ATA-CB-CB2RA	Anti AP 1	0.1 ML	Human Rat	Rabbit		IHC(f) IHC(p) WB IF LM	polyclonal		Serum
ATA-CB-CB3RA	Anti AP 2	0.1 ML	Human Rat	Rabbit		IHC(f) IHC(p) WB IF LM	polyclonal		Serum
PPX-PP-H7507-00	Anti AR	0.1 ML	Human Rat Mouse	Mouse	H7507	WB ELISA IHC(p) IP	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
CAC-SDT-01-APX	Anti Ascorbate Peroxidase	200 UL	Plant	Rabbit		WB IHC	polyclonal	IgG	Ig fraction - Ammonium Sulphate
CAC-TMD-PH-AT12	Anti Atg12	100 UL	Mouse	Rabbit		WB IHC	polyclonal	IgG	Serum
CAC-TMD-PH-AT16L	Anti Atg16L	100 UL	Mouse Rat Human	Rabbit		WB IHC	polyclonal	IgG	Serum
YII-YP060-EX	Anti ATP Receptor:P2X4 (370-388)	50 UL	Rat	Rabbit		IHC	polyclonal		Serum
KAL-KR051	Anti Aurora A	25 UG	Human	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
KAL-KM105	Anti Autotaxin /ENPP2	25 UG	-	Rabbit		IHC	polyclonal		Affinity Purified
KAL-KW381	Anti Avidin HRP	1 MG		Rabbit		WB ELISA IHC IC	polyclonal		
KAL-KE029	Anti b (0,+)-Type Amino Acid Transporter 1 (BAT1)	200 UG	Rat	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
KAL-KR093	Anti BNAS2	25 UG	Mouse	Rabbit		WB IHC IP	polyclonal	IgG	Affinity Purified

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
YII-Y170-EX	Anti Bombesin	50 UL	Porcine Canine Monkey Guinea Pig	Guinea Pig		IHC RIA	polyclonal		Serum
LSL-LB-4335	Anti Bone Sialoprotein	100 UL	Human Mouse Rat Bovine	Rabbit		WB ELISA IHC	polyclonal		Serum
YII-YP050-EX	Anti Brain NOS-1 (998-1024)	50 UL	Human	Rabbit		IHC	polyclonal		Serum
HKD-AB43	Anti branched-chain $\alpha$ -keto acid dehydrogenase (BCKDH) E2	150 UL	Rat	Mouse		WB IHC IP	monoclonal	IgG1	
HKD-AB44	Anti branched-chain $\alpha$ -keto acid dehydrogenase (BCKDH) kinase	100 UL	Rat	Mouse		WB IHC IP	monoclonal	IgG1	
YII-Y221-EX	Anti C Peptide I	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
YII-Y220-EX	Anti C Peptide II	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
KAL-KR084	Anti Calpain 6	25 UG	Human Mouse	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
KAL-KW157	Anti Calretinin	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p)	polyclonal	IgG	Affinity Purified
KAL-KY041	Anti CaM Kinase II ( $\delta$ 1-4)	200 UG	Human Mouse Rat	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
YII-Y470-EX	Anti Carassius RFamide	50 UL		Rabbit		IHC(p) IHC(f) ELISA			
KAL-KR066	Anti Cation Chloride Cotransporter 9 (CCC9)	25 UG	Mouse	Rabbit		IHC	polyclonal	IgG	Affinity Purified
YII-YP030-EX	Anti CCK-8	50 UL	Human Rat	Rabbit		IHC RIA	polyclonal		Serum
SID-IDL011	Anti CD20	100 UG	Human	Mouse	L26	WB IHC	monoclonal	IgG2a	

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
SID-IDL014	Anti CD20 Biotin	50 UG	Human	Mouse	L26	WB IHC	monoclonal	IgG2a	
KAL-KS131	Anti CD298	50 UG	Human	Mouse	5G10	FC IP IHC	monoclonal	IgG2ak	Ig fraction - Protein G
KAL-KN141	Anti CD3	50 UG		Mouse	NU-T3	IHC FC	monoclonal	IgG2a	Ig fraction - Protein G
KAL-KW162	Anti CD40L	100 UG	Human Mouse Rat	Rabbit		WB IHC(p)	polyclonal	IgG	Affinity Purified
KAL-KS124	Anti CD9	1 VIAL	Human	Mouse	6D11	WB IHC FC IP	monoclonal	IgG2ak	Ig fraction - Protein G
KAL-KH025	Anti CEL	50 UG	-	Mouse	KNH-30	IHC ELISA	monoclonal	IgG1	Affinity Purified
KAL-KH025-01	Anti CEL Biotin	50 UG	-	Mouse	KNH-30	IHC ELISA	monoclonal	IgG1	Affinity Purified
KAL-KH025-02	Anti CEL HRP	50 UG	-	Mouse	KNH-30	IHC ELISA	monoclonal	IgG1	Affinity Purified
ALF-274042003	Anti Ceramide	0.25 MG		Mouse	NHCER-2	IHC IB	monoclonal	IgM	Ig fraction - Ion Exchange /Gel Filtration
YII-Y340-EX	Anti CGRP	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
KAL-KR072	Anti Chemokine-like Factor 1(CKLF1)	25 UG	Human	Rabbit		IHC	polyclonal	IgG	Affinity Purified
YII-Y292-EX	Anti Chromogranin A (343-360)	50 UL	Rat	Rabbit		IHC(f) RIA	polyclonal		Serum
YII-Y293-EX	Anti Chromogranin A (344-374)	50 UL	Human	Rabbit		IHC(f) EIA RIA	polyclonal		Serum
YII-Y290-EX	Anti Chromogranin A (359-389)	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y291-EX	Anti Chromogranin A (94-130)	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
ATA-CB-CB1RA	Anti Clathrin	0.1 ML	Rat	Rabbit		IHC(f) IHC(p) IF LM	polyclonal		Serum
KAL-KH011	Anti CML	50 UG	-	Mouse	CMS-10	IHC	monoclonal	IgG1	Affinity Purified

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
KAL-KH024	Anti CML	50 UG	-	Mouse	NF-1G	IHC ELISA	monoclonal	IgG2a	Affinity Purified
KAL-KH024-01	Anti CML	50 UG	-	Mouse	NF-1G	IHC	monoclonal	IgG2a	Affinity Purified
KAL-KH011-01	Anti CML Biotin	50 UG	-	Mouse	CMS-10	IHC	monoclonal	IgG1	Affinity Purified
KAL-KH011-02	Anti CML HRP	50 UG	-	Mouse	CMS-10	IHC ELISA	monoclonal	IgG1	Affinity Purified
KAL-KH024-02	Anti CML HRP	50 UG	-	Mouse	NF-1G	IHC	monoclonal	IgG2a	Affinity Purified
LSL-LB-1102	Anti Collagen 1	100 UL	Human Mouse Rat Bovine Porcine Guinea Pig Mink	Rabbit		WB IHC	polyclonal		Serum
LSL-LB-1190	Anti Collagen 1	100 UL	Human Bovine Porcine Monkey Guinea Pig	Rabbit		IHC	polyclonal		Serum
LSL-LB-1196	Anti Collagen 1	100 UL	Human Bovine Porcine Canine	Rabbit		IHC	polyclonal		Serum
LSL-LB-1197	Anti Collagen 1	100 UL	Human Bovine Porcine Sheep Guinea Pig	Rabbit		IHC	polyclonal		Serum
LSL-LB-0092	Anti Collagen 10	100 UL	Human Mouse Rat Chicken	Rabbit		IHC	polyclonal		Serum
LSL-LB-1297	Anti Collagen 2	100 UL	Human Mouse Rat Bovine Guinea Pig	Rabbit		IHC	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
LSL-LB-1300	Anti Collagen 3	100 UL	Human Bovine Porcine Mink	Rabbit		IHC	polyclonal		Serum
LSL-LB-1387	Anti Collagen 3	100 UL	Human Mouse Rat Bovine Guinea Pig	Rabbit		IHC	polyclonal		Serum
LSL-LB-1393	Anti Collagen 3	100 UL	Human Mouse Rat Bovine Porcine	Rabbit		IHC	polyclonal		Serum
DFK-F-59-EX	Anti Collagen 4	1 ML	Human	Mouse	IV-4H12	WB IHC(p) IHC(f) EIA	monoclonal	IgG1κ	Purified
LSL-LB-0445	Anti Collagen 4	100 UL	Human Mouse Rat Bovine Monkey	Rabbit		WB IHC	polyclonal		Serum
LSL-LB-1403	Anti Collagen 4	100 UL	Mouse Rat Guinea Pig Mink	Rabbit		IHC	polyclonal		Serum
LSL-LB-1407	Anti Collagen 4	100 UL	Human Mouse Rat Bovine Porcine Goat Guinea Pig Mink	Rabbit		IHC	polyclonal		Serum
LSL-LB-1503	Anti Collagen 5	100 UL	Human Mouse Rat Bovine	Rabbit		IHC	polyclonal		Serum
LSL-LB-1597	Anti Collagen 5	100 UL	Human Bovine Porcine Goat Guinea Pig	Rabbit		IHC	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
LSL-LB-1697	Anti Collagen 6	100 UL	Human Mouse Rat Bovine Guinea Pig	Rabbit		IHC	polyclonal		Serum
LSL-LB-0771	Anti Collagen 7	100 UL	Human Mouse Rat Bovine	Rabbit		IHC	polyclonal		Serum
LSL-LB-0883	Anti Collagen 8	100 UL	Human Mouse Rat Bovine	Rabbit		IHC	polyclonal		Serum
LSL-LB-0903	Anti Collagen 9	100 UL	Human Mouse Rat Bovine	Rabbit		IHC	polyclonal		Serum
KAL-KW167	Anti Connexin43	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p) IC	polyclonal	IgG	Affinity Purified
PPX-PP-H8132-00	Anti COUP-TF 1	0.1 ML	Human Mouse Rat	Mouse	H8132	WB ELISA IHC IP	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-H7147-00	Anti COUP-TF 2	0.1 ML	Human Mouse Rat	Mouse	H7147	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-H8124-00	Anti COUP-TF I	0.1 ML	Human Mouse Rat	Mouse	H8124	ELISA WB IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
CAC-NM-DND-001	Anti CPDs	1 VIAL		Mouse	TDM-2	ELISA IC IHC	monoclonal	IgG2ak	Semi Purified
YII-Y222-EX	Anti C-Peptide I	50 UL	Mouse	Rabbit		IHC(p) IHC(f) IA	polyclonal		Serum
YII-Y223-EX	Anti C-Peptide II	50 UL	Mouse	Rabbit		IHC(p) IHC(f) IA	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
YII-Y210-EX	Anti CRF (24-41)	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
YII-Y211-EX	Anti CRF (3-41)	50 UL	Human Mouse Rat	Rabbit		IHC(f) RIA	polyclonal		Serum
NOF-N213630-EX	Anti Crotonaldehyde	30 UG	-	Mouse		IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
KAL-KW169	Anti CRP	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p)	polyclonal	IgG	Affinity Purified
LSL-LB-5201	Anti Cry j-1	100 UL		Rabbit		WB IHC(p) IHC(f) ELISA	polyclonal		Serum
LSL-LB-5202	Anti Cry j-2	100 UL		Rabbit		WB IHC(p) IHC(f) ELISA	polyclonal		Serum
CAC-PUH-011N-A	Anti Cu / Zn-SOD	100 UL	Human		5C10	WB IHC	monoclonal	IgG1k	Supernatant
CAC-PUH-011N-B	Anti Cu / Zn-SOD	500 UL	Human		5C10	WB IHC	monoclonal	IgG1k	Supernatant
KAL-KE021	Anti Cystine / Glutamic Acid Transporter (xCT)	25 UG	Mouse	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
CAC-TNL-001-CAM	Anti dCAMKII	100 UL	Drosophila	Mouse	18	WB IHC	monoclonal	IgG2b	Ascites
KAL-KR104	Anti DC-STAMP(FIND)	25 UG	Mouse	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KG410	Anti Dcun1d4	25 UG	Mouse	Rabbit		IHC(p)	polyclonal		Affinity Purified
NNS-MBY-020P-EX	Anti Dibromo-Tyrosine (DiBrY)	20 UG		Mouse	3A5	IHC WB	monoclonal	IgG1k	Ig fraction - Protein A
NNS-MDT-020P-EX	Anti Dityrosine (DT)	20 UG		Mouse	1C3	IHC WB	monoclonal	IgG2ak	Ig fraction - Protein A
PPX-PP-H9929A-00	Anti EAR2	0.1 ML	Human Rat Mouse	Mouse	H9929A	WB ELISA IHC	monoclonal		
PPX-PP-N2025-00	Anti EAR2	0.1 ML	Human Rat Mouse	Mouse	N2025	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
KAL-KW172	Anti EDNR-B	100 UG	Human Mouse Rat	Rabbit		WB IHC(p) IC	polyclonal	IgG	Affinity Purified
YII-Y230-EX	Anti EGF	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
YII-Y231-EX	Anti EGF	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
LSL-LB-2187	Anti Elastin	100 UL	Human Mouse Rat Bovine	Rabbit		WB IHC	polyclonal		Serum
PPX-PP-H5844-00	Anti ERR $\alpha$	0.1 ML	Human Rat Mouse	Mouse	H5844	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-H6705-00	Anti ERR $\beta$	0.1 ML	Human Rat Mouse	Mouse	H6705	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-H6812-00	Anti ERR $\gamma$	0.1 ML	Human Rat Mouse	Mouse	H6812	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-H4624-00	Anti ER $\alpha$	0.1 ML	Human Rat Mouse	Mouse	H4624	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
CAC-OUA-P01-A	Anti EVC	100 UL	Rat Mouse	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
SML-ROA003-EX	Anti F3	200 UL	Mouse	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
KYO-98140	Anti Fc $\epsilon$ R1	200 UG	Human	Mouse	CRA1	WB IHC(p) IHC(f) FC	monoclonal	IgG2b	
BAM-72-001	Anti Fc $\epsilon$ R1 $\alpha$	100 UG	Human	Mouse	CRA1	WB IHC FC	monoclonal	IgG2bk	Purified
BAM-72-005	Anti Fc $\epsilon$ R1 $\alpha$	100 UG	Human	Mouse	CRA2	WB IHC FC	monoclonal	IgG1 $\kappa$	Purified
BAM-72-003	Anti Fc $\epsilon$ R1 $\alpha$ Biotin	50 UG	Human	Mouse	CRA1	WB FC IHC	monoclonal	IgG2bk	
BAM-72-007	Anti Fc $\epsilon$ R1 $\alpha$ Biotin	50 UG	Human	Mouse	CRA2	WB FC IHC	monoclonal	IgG1 $\kappa$	
BAM-72-004	Anti Fc $\epsilon$ R1 $\alpha$ FITC	50 UG	Human	Mouse	CRA1	WB FC IHC	monoclonal	IgG2bk	
BAM-72-008	Anti Fc $\epsilon$ R1 $\alpha$ FITC	50 UG	Human	Mouse	CRA2	FC IHC	monoclonal	IgG1 $\kappa$	
KAL-KX411	Anti FFAR2	50 UG	Human Rat	Mouse	20D2	IHC ELISA	monoclonal	IgG2bk	Ig fraction - Protein G
YII-Y250-EX	Anti FGF Basic	50 UL	Human Bovine	Rabbit		IHC(f) RIA	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
YII-Y251-EX	Anti FGF Basic (1-9)	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
LSL-LB-2297	Anti Fibrillin	100 UL	Human Mouse Rat Bovine	Rabbit		WB IHC	polyclonal		Serum
LSL-LB-1027	Anti Fibronectin	100 UL	Human Mouse Rat Bovine Goat Guinea Pig	Rabbit		IHC	polyclonal		Serum
PPX-PP-A9033A-00	Anti FXR	0.1 ML	Human Mouse Rat	Mouse	A9033A	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
YII-YP071-EX	Anti G Prptein:Gα i3 (105-127)	50 UL	Human	Rabbit		WB IHC	polyclonal		Serum
YII-YP070-EX	Anti G Prptein:Gα oIF (23-38)	50 UL	Human	Rabbit		WB IHC	polyclonal		Serum
PPX-PP-A8620A-00	Anti G9a(H3-K9-HMTase)	0.1 ML	Mouse Human	Mouse	A8620A	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
KAL-KW177	Anti GAD67	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p)	polyclonal	IgG	Affinity Purified
YII-Y186-EX	Anti Galanin	50 UL	Rat	Rabbit		WB IHC	polyclonal		Serum
YII-Y182-EX	Anti Galanin (1-15)	50 UL	Rat Tuna	Rabbit		IHC RIA	polyclonal		Serum
YII-Y183-EX	Anti Galanin (15-29)	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y184-EX	Anti Galanin (9-30)	50 UL	Human Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y312-EX	Anti GAP (28-56)	50 UL	Rat	Rabbit		IHC(f) RIA	polyclonal		Serum
YII-Y311-EX	Anti GAP (34-56)	50 UL	Human Rat Porcine	Rabbit		IHC(f) RIA	polyclonal		Serum
YII-Y111-EX	Anti Gastrin	50 UL	Rat Canine	Rabbit		IHC RIA	polyclonal		Serum
YII-Y110-EX	Anti Gastrin34 (1-15)	50 UL	Human Rat	Rabbit		IHC RIA	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
HKD-AB36	Anti Gelatinase B	100 UL	Killifish	Rabbit		WB IHC(f) IP	polyclonal	-	-
BAM-73-003	Anti GENA (clone TRA98)	50 UG		Mouse	TRA98	WB IHC			
BAM-73-004	Anti GENA (clone TRA98)	250 UG		Mouse	TRA98	WB IHC			
SML-ROI003-EX	Anti GFAP	0.2 ML	Human Mouse Rat	Rabbit		WB IHC(p)	polyclonal	IgG	Purified
KAL-KR069	Anti Ghrelin	25 UG	Rat	Rabbit		IHC(p)	polyclonal	IgG	Affinity Purified
YII-Y390-EX	Anti GHRH	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y101-EX	Anti GIP	50 UL	Human Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y102-EX	Anti GIP	50 UL	Rat	Rabbit		IHC(p) IHC(f) EIA	polyclonal		Serum
YII-Y100-EX	Anti GIP (1-30)·OH	50 UL	Rat Porcine	Rabbit		IHC RIA	polyclonal		Serum
YII-Y103-EX	Anti GIP (18-42)	50 UL	Rat	Rabbit		IHC(p) IHC(f)	polyclonal		Serum
YII-Y324-EX	Anti Glicentin (1-32)	50 UL	Rat	Rabbit		IHC(p) IHC(f) EIA IA	polyclonal		Serum
YII-Y320-EX	Anti GLP 1 (7-36)amide	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
YII-Y322-EX	Anti GLP 2	50 UL	Rat	Rabbit		IHC(p) EIA	polyclonal		Serum
YII-Y321-EX	Anti GLP 2 (14-33)	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
PPX-PP-B0422-00	Anti GLP(Eu-HMTase)	0.1 ML	Mouse Human	Mouse	B0422	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
YII-Y323-EX	Anti GLP-2	50 UL	Mouse Rat	Rabbit		IHC(p) IHC(f) EIA IA	polyclonal		Serum
ATA-CB-GP1RA	Anti Glutathione peroxidase (GPx 1)	0.1 ML	Human Rat	Rabbit		IHC(f) IHC(p) WB ELISA IF LM	polyclonal		Serum
ATA-CB-GP1SA	Anti Glutathione peroxidase (GPx 1)	0.1 ML	Human Rat	Sheep		IHC(f) IHC(p) WB ELISA IF LM	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
ATA-CB-GPPRA	Anti Glutathione peroxidase (GPx P)	0.1 ML	Human Rat	Rabbit		IHC(f) IHC(p) WB ELISA IF LM	polyclonal		Serum
KAL-KX412	Anti GPR39	50 UG	Human Rat	Mouse	8H11	IHC ELISA	monoclonal	IgG2ak	Ig fraction - Protein G
PPX-PP-H8004-00	Anti GR common	0.1 ML	Human Rat	Mouse	H8004	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-H8031-00	Anti GR common	0.1 ML	Human Rat Mouse	Mouse	H8031	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
KAL-KW181	Anti GST $\pi$	100 UG	Human Mouse Rat	Rabbit		WB IHC(p) IC	polyclonal	IgG	Affinity Purified
KAL-KG408	Anti Haa0	25 UG	Mouse	Rabbit		IHC(p)	polyclonal		Affinity Purified
SIM-2ZHB16	Anti HBs determinant a	0.5 ML		Mouse	5124A	ELISA WB IHC	monoclonal	IgM	Purified
HKD-AB-46	Anti HLA- A, B, C	100 UL	Human	Mouse	EMR8-5	WB IHC	monoclonal	IgG1	-
SID-IDL021	Anti HLA A,B,C	100 UG	Human	Mouse	EMR8-5	WB IHC	monoclonal	IgG1	
SID-IDL024	Anti HLA A,B,C Biotin	50 UG	Human	Mouse	EMR8-5	WB IHC	monoclonal	IgG1	
HKD-AB-46-B	Anti HLA class I -A,B,C Biotin	100 UL	Human	Mouse	EMR8-5	IHC WB	monoclonal	IgG1	
HKD-AB-46-F	Anti HLA class I -A,B,C FITC	100 UL	Human	Mouse	EMR8-5	IHC WB	monoclonal	IgG1	
HKD-AB-46-H	Anti HLA class I -A,B,C HRP	100 UL	Human	Mouse	EMR8-5	WB IHC	monoclonal	IgG1	
KAL-KR077	Anti Hlxb9	25 UG	Mouse	Rabbit		IHC	polyclonal	IgG	Affinity Purified
PPX-PP-N3224-00	Anti HNF4 $\gamma$	0.1 ML	Human Mouse Rat	Mouse	N3224	WB ELISA IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-K9218-00	Anti HNF4 $\alpha$ (isoforms 1-6)	0.1 ML	Human Mouse Rat	Mouse	K9218	WB ELISA IP Gel Shift IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-H1415-00	Anti HNF4 $\alpha$ (isoforms 1,2,4,5,7 and 8)	0.1 ML	Human Mouse Rat	Mouse	H1415	WB ELISA IP Gel Shift IHC ChIP	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-H6939-00	Anti HNF4 $\alpha$ 7(isoforms 7-9)	0.1 ML	Human Rat Mouse	Mouse	H6939	WB ELISA IP IHC	monoclonal	IgG1	Ig fraction - Ammonium Sulphate

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
PPX-PP-B6502A-00	Anti HNF4 $\gamma$	0.1 ML	Human Rat	Mouse	B6502A	ELISA WB IP IHC	monoclonal	IgG1	Ig fraction - Ammonium Sulphate
CAC-RIK-B-HP	Anti Homer (Pan)	50 UL	Human Mouse Rat	Rabbit		WB IHC	polyclonal	IgG	
CAC-RIK-B-H1	Anti Homer1	50 UL	Human Mouse Rat	Rabbit		WB IHC	polyclonal	IgG	
CAC-RIK-B-H2	Anti Homer2	50 UL	Human Mouse Rat	Rabbit		WB IHC	polyclonal	IgG	
CAC-RIK-B-H3	Anti Homer3	50 UL	Human Mouse Rat	Rabbit		WB IHC	polyclonal	IgG	
SIM-2ZHPUA	Anti HP Urease A	0.5 ML	Helicobacter pylori	Mouse	K2127B	WB IHC(f) ELISA	monoclonal	IgG1	Purified
SIM-2ZHPUB	Anti HP Urease B	0.5 ML	Helicobacter pylori	Mouse	K2031B	WB IHC(f) ELISA	monoclonal	IgG1	Purified
BAM-70-221	Anti HP1 $\alpha$ /CBX5	50 UG	Human	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
BAM-70-223	Anti HP1 $\beta$ /CBX1	50 UG	Human Hamster	Rabbit		WB IHC ChIP	polyclonal	IgG	Affinity Purified
BAM-70-225	Anti HP1 $\gamma$ /CBX3	50 UG	Human Hamster	Rabbit		WB IHC ChIP	polyclonal	IgG	Affinity Purified
KAL-KW240	Anti HSP27	100 UG	Human Rat	Rabbit		WB IHC(p)	polyclonal	IgG	Affinity Purified
KAL-KW228	Anti HSP60	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p) IHC(f) IC	polyclonal	IgG	Affinity Purified
YII-Y440-EX	Anti Humanin	50 UL		Rabbit		IHC EIA	polyclonal		Serum
KAL-KW372	Anti IgA Biotin	1 MG	Human	Goat		WB ELISA IHC IC	polyclonal	IgG	Affinity Purified
KAL-KW394	Anti IgA FITC	1 MG	Human	Goat		IHC FC IC IF	polyclonal	IgG	Affinity Purified
KAL-KW378	Anti IgA HRP	1 MG	Human	Goat		WB ELISA IHC IC	polyclonal	IgG	Affinity Purified

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
CAC-YCU-MK-TF01	Anti IGF1P-rP1	100 UG	Human	Mouse	88	WB ELISA IHC	monoclonal	IgG	Affinity Purified
KAL-KW390	Anti IgG FITC	1 MG	Mouse	Goat		IHC FC IC IF	polyclonal	IgG	Affinity Purified
KAL-KW391	Anti IgG FITC	1 MG	Rabbit	Goat		IHC FC IC IF	polyclonal	IgG	Affinity Purified
KAL-KW392	Anti IgG FITC	1 MG	Rat	Rabbit		IHC FC IC IF	polyclonal	IgG	Affinity Purified
KAL-KW393	Anti IgG FITC	1 MG	Goat	Rabbit		IHC FC IC IF	polyclonal	IgG	Affinity Purified
KAL-KW395	Anti IgG FITC	1 MG	Human	Rabbit		IHC FC IC IF	polyclonal	IgG	Affinity Purified
KAL-KW382	Anti IgG GCP	1 MG	Mouse	Goat		IHC	polyclonal	IgG	Affinity Purified
KAL-KW383	Anti IgG GCP	1 MG	Mouse	Goat		IHC	polyclonal	IgG	Affinity Purified
KAL-KW384	Anti IgG GCP	1 MG	Rabbit	Goat		IHC	polyclonal	IgG	Affinity Purified
KAL-KW385	Anti IgG GCP	1 MG	Rabbit	Goat		IHC	polyclonal	IgG	Affinity Purified
KAL-KW386	Anti IgG GCP	1 MG	Goat	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KW387	Anti IgG GCP	1 MG	Goat	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KW388	Anti IgG GCP	1 MG	Human	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KW389	Anti IgG GCP	1 MG	Human	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KW374	Anti IgG HRP	1 MG	Mouse	Goat		WB ELISA IHC IC	polyclonal	IgG	Affinity Purified
KAL-KW375	Anti IgG HRP	1 MG	Rabbit	Goat		WB ELISA IHC IC	polyclonal	IgG	Affinity Purified
KAL-KW376	Anti IgG HRP	1 MG	Rat	Rabbit		WB ELISA IHC IC	polyclonal	IgG	Affinity Purified
KAL-KW377	Anti IgG HRP	1 MG	Goat	Rabbit		WB ELISA IHC IC	polyclonal	IgG	Affinity Purified
KAL-KW379	Anti IgG HRP	1 MG	Human	Rabbit		WB ELISA IHC IC	polyclonal	IgG	Affinity Purified

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
KAL-KW397	Anti IgG TRITC	1 MG	Mouse	Goat		IHC IC IF	polyclonal	IgG	Affinity Purified
KAL-KW398	Anti IgG TRITC	1 MG	Rabbit	Goat		IHC IC IF	polyclonal	IgG	Affinity Purified
KAL-KW399	Anti IgG TRITC	1 MG	Goat	Rabbit		IHC IC IF	polyclonal	IgG	Affinity Purified
KAL-KW400	Anti IgG TRITC	1 MG	Human	Rabbit		IHC IC IF	polyclonal	IgG	Affinity Purified
KAL-KW373	Anti IgM Biotin	1 MG	Human	Goat		WB ELISA IHC IC	polyclonal	IgG	Affinity Purified
KAL-KW396	Anti IgM FITC	1 MG	Human	Goat		IHC FC IC IF	polyclonal	IgG	Affinity Purified
KAL-KW380	Anti IgM HRP	1 MG	Human	Goat		WB ELISA IHC IC	polyclonal	IgG	Affinity Purified
KAL-KR070	Anti Initiate Factor 3 (IF3)	25 UG	Mouse	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KR079	Anti Initiation Factor 5 (IF5)	25 UG	Mouse	Rabbit		IHC	polyclonal	IgG	Affinity Purified
YII-Y370-EX	Anti Insulin	50 UL	Human	Guinea Pig		IHC RIA	polyclonal		Serum
KAL-KW186	Anti ITGA1	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p)	polyclonal	IgG	Affinity Purified
KAL-KW187	Anti ITGA4	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p) IC	polyclonal	IgG	Affinity Purified
LSL-LB-1003	Anti Keratin	100 UL	Human Mouse Rat Bovine Sheep Mink Sea Lion	Rabbit		IHC	polyclonal		Serum
KAL-KR096	Anti KF-1	25 UG	Mouse	Rabbit		IHC	polyclonal	IgG	Affinity Purified
LSL-LB-1013	Anti Laminin	100 UL	Human Mouse Rat Bovine Porcine Guinea Pig	Rabbit		IHC	polyclonal		Serum
KAL-KG405	Anti Leprotl1	25 UG	Mouse	Rabbit		IHC(p)	polyclonal		Affinity Purified
YII-YC040-EX	Anti Leptin	50 UL	Rat	Rabbit		IHC(p) IHC(f) RIA	polyclonal		Serum
YII-YP010-EX	Anti LH-RH	50 UL	Human Rat	Rabbit		IHC RIA	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
KAL-KE023	Anti L-Type Amino Acid Transporter 1 (LAT1)	20 UG	Human	Mouse	4D9	IHC WB	monoclonal	IgM	Affinity Purified
KAL-KE026	Anti L-Type Amino Acid Transporter 1 (LAT1)	25 UG	Human	Rabbit		IHC WB	polyclonal	IgG	Affinity Purified
YII-Y380-EX	Anti Luminal CCK-RF	50 UL	Rat	Rabbit		IHC(p) IHC(f) RIA	polyclonal		Serum
PPX-PP-PPZ0412-00	Anti LXR $\alpha$ Ligand Binding Domain	0.1 ML	Human Rat Mouse	Mouse	PPZ0412	WB ELISA IP Gel Shift IHC ChIP	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
KAL-KT014	Anti Macrophage / Dendritic Cells	50 UG	Rat	Mouse	RM-4	IHC	monoclonal	IgG1	Affinity Purified
KAL-KT015	Anti Macrophage / Monocyte	50 UG	Rabbit	Mouse	RbM-2	IHC	monoclonal	IgG1	Affinity Purified
KAL-KT118	Anti Macrophage Scavenger Receptor A	50 UG	Human	Mouse	SRA-C6	WB IHC	monoclonal	IgG1	Affinity Purified
KAL-KT118-S	Anti Macrophage Scavenger Receptor A (additive-free)	50 UG	Human	Mouse	SRA-C6	WB IHC	monoclonal	IgG1	Affinity Purified
KAL-KT013	Anti Macrophage Surface Antigen	50 UG	Human	Mouse	AM-3K	IHC	monoclonal	IgG1	Affinity Purified
NOF-N213530-EX	Anti Malondialdehyde	30 UG	-	Mouse	1F83	IHC	monoclonal	IgG2a $\lambda$	Ig fraction - Protein A
KAL-KR083	Anti Matrix Gla Protein (MGP)	50 UG	Mouse	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KW191	Anti MBP	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p) IC	polyclonal	IgG	Affinity Purified
SML-ROA001-EX	Anti MBP, phospho	200 UL	Human Mouse Rat	Rabbit		WB IHC(p) ELISA DB	polyclonal	IgG	Ig fraction - Protein A
KAL-KA009	Anti Metallothionein	100 UG	Human Mouse Rat Rabbit	Mouse	1A12	WB IHC	monoclonal	IgG1	Affinity Purified
KAL-KA009-01	Anti Metallothionein Biotin	100 UG	Human Mouse Rat Rabbit	Mouse	1A12	WB IHC	monoclonal	IgG1	Affinity Purified

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
YII-Y140-EX	Anti Met-Enkephalin-Arg-Gly-Leu	50 UL	Human Rat Bovine Canine Monkey Guinea Pig Feline	Rabbit		IHC RIA	polyclonal		Serum
NOF-N213430-EX	Anti Methylglyoxal	30 UG	-	Mouse		IHC	monoclonal	IgG2a	Ig fraction - Protein A
KAL-KJ091	Anti MFH	50 UG	Rat	Mouse	A3	WB IHC	monoclonal	IgG1	Affinity Purified
DFK-F-67-EX	Anti MMP 1	1 ML	Human Rabbit	Mouse	41-1E5	WB IHC(p) IHC(f)	monoclonal	IgG2ak	Purified
DFK-F-89-EX	Anti MMP 13	1 ML	Human	Mouse	181-15A12	WB IHC(p) EIA	monoclonal	IgG1k	Purified
DFK-F-94-EX	Anti MMP 13	1 ML	Human	Mouse	181-14G11	WB IHC(p)	monoclonal	IgG1k	
DFK-F-68-EX	Anti MMP 2	1 ML	Human Mouse Rat Bovine Rabbit	Mouse	42-5D11	WB IHC(p)	monoclonal	IgG1k	Purified
DFK-F-73-EX	Anti MMP 2	1 ML	Human Mouse Bovine	Mouse	75-7F7	WB IHC(p) IHC(f) EIA	monoclonal	IgG1k	Purified
DFK-F-66-EX	Anti MMP 3	1 ML	Human Rat Porcine Sheep	Mouse	55-2A4	WB IHC(p) IHC(f) EIA	monoclonal	IgG1k	Purified
DFK-F-83-EX	Anti MMP 8	1 ML	Human	Mouse	115-13D2	WB IHC(p)	monoclonal	IgG1k	Purified
DFK-F-69-EX	Anti MMP 9	1 ML	Human Rat Rabbit Guinea Pig	Mouse	56-2A4	WB IHC(p) IHC(f)	monoclonal	IgG1k	Purified
CAC-PUH-021N-A	Anti MnSOD	100 UL	Human		2C6	WB IHC	monoclonal	IgM	Supernatant

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
CAC-PUH-021N-B	Anti MnSOD	500 UL	Human		2C6	WB IHC	monoclonal	IgM	Supernatant
KAL-KR085	Anti MOK	25 UG	Mouse	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
YII-Y120-EX	Anti Motilin	50 UL	Canine	Rabbit		IHC RIA	polyclonal		Serum
PPX-PP-H3122-00	Anti MR	0.1 ML	Human Mouse Rat	Mouse	H3122	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
KAL-KT022	Anti MSR-A (CD204)	50 UG	Human Mouse Rat Bovine Rabbit Porcine Canine Monkey	Mouse	SRA-E5	WB IHC	monoclonal	IgG1	Affinity Purified
LSL-LB-3004	Anti Myelin protein zero (PO)	100 UL	Mouse Rat Chicken	Rabbit		WB IHC IF	polyclonal		
KAL-KG409	Anti Mylpf	25 UG	Mouse	Rabbit		IHC(p)	polyclonal		Affinity Purified
LSL-LB-1074	Anti Myosin	100 UL	Chicken	Rabbit		IHC	polyclonal		Serum
YII-Y400-EX	Anti Na Channel 2 (467-485)	50 UL	Rat	Rabbit		IHC(f)	polyclonal		Serum
LSL-LB-6198	Anti Na/ K ATPase	100 UL	Mouse Rat Bovine Guinea Pig Equine	Rabbit		IHC(p) IHC(f)	polyclonal		Serum
IKG-KMSLS3-1	Anti Native $\alpha$ 1 Antitrypsin	100 UG	Human	Mouse	S3-1	EIA WB IHC(f)	monoclonal	IgG1	
BAM-73-100	Anti Nestin	200 UG	Mouse	Rat	7A3	IHC IC	monoclonal	IgG2bk	Purified
CAC-GU01-M02AS-A	Anti Neurofilament H/M ,phospho	50 UL		Mouse	NFP-1D	IHC	monoclonal	IgM	Ascites

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
YII-Y130-EX	Anti Neurotensin	50 UL	Human Rat Bovine Monkey	Rabbit		IHC RIA	polyclonal		Serum
PPX-PP-N1404-00	Anti NGFI -B $\beta$	0.1 ML	Human Rat Mouse	Mouse	N1404	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
KAL-KW199	Anti NMDAR-2A	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p) IC	polyclonal	IgG	Affinity Purified
KAL-KW200	Anti NMDAR-2B	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p) IC	polyclonal	IgG	Affinity Purified
YII-YP051-EX	Anti NOS 2(77-105)	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y060-EX	Anti NPY	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
YII-Y061-EX	Anti NPY (1-19)	50 UL	Human Mouse Rat	Rabbit		IHC	polyclonal		Serum
KAL-KW203	Anti NT-3	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p)	polyclonal	IgG	Affinity Purified
BAM-70-310	Anti Nup98	200 UG	Human Mouse Rat		2H10	WB IHC ELISA DB	monoclonal	IgG2ck	Purified
YII-Y460-EX	Anti Obestatin	50 UL	Mouse Rat	Rabbit		IHC EIA	polyclonal		Serum
YII-Y461EX	Anti Obestatin	50 UL	Human	Rabbit		IHC	polyclonal		Serum
CAC-RIK-B-OP	Anti Opalin	100 UG	Mouse Rat	Rabbit		WB IHC IP	polyclonal	IgG	
YII-Y450-EX	Anti Orexin A	50 UL	Human Mouse Rat Bovine	Rabbit		IHC ELISA	polyclonal	IgG	Serum
YII-Y451-EX	Anti Orexin B	50 UL	Mouse Rat	Rabbit		IHC ELISA	polyclonal	IgG	Serum
KAL-KE031	Anti Organic Anion Transporter 2 (OAT2)	25 UG	Human	Rabbit		IHC	polyclonal	IgG	Affinity Purified

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
KAL-KE032	Anti Organic Anion Transporter 3 (OAT3)	25 UG	Human	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KE035	Anti Organic Anion Transporter 3 (OAT3)	25 UG	Rat	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KE033	Anti Organic Anion Transporter 4 (OAT4)	20 UG	Human	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KE038	Anti Organic Anion Transporter1	25 UG	Human	Rabbit		IHC	polyclonal		Affinity Purified
KAL-KR071	Anti Organic Cation Transporter 3 (OCT3)	25 UG	Rat	Rabbit		IHC	polyclonal	IgG	Affinity Purified
YII-Y430-EX	Anti Orphanin FQ / Nociceptin	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
LSL-LB-4005	Anti Osteocalcin	100 UL	Human Mouse Rat Bovine	Rabbit		WB IHC(p) IHC(f)	polyclonal		Serum
LSL-LB-4115	Anti Osteonectin	100 UL	Human Mouse Rat Bovine	Rabbit		WB IHC(p) IHC(f)	polyclonal		Serum
LSL-LB-4225	Anti Osteopontin	100 UL	Human Mouse Rat Bovine	Rabbit		WB IHC(p) IHC(f)	polyclonal		Serum
IKG-KMSL3F4	Anti Oxidized $\alpha$ 1 Antitrypsin	100 UG	Human	Mouse	3F4	EIA WB IHC(f)	monoclonal	IgG1	
KAL-KX413	Anti P2RY12	50 UG	Human Rat	Mouse	2A6	IHC ELISA	monoclonal	IgG2bk	Ig fraction - Protein G
KAL-KW205	Anti P73 $\alpha$	100 UG	Human Mouse Rat	Rabbit		WB IHC(p)	polyclonal	IgG	Affinity Purified
ALF-274042065	Anti PA	0.20 MG		Mouse	KA-1	ELISA IHC IB	monoclonal	IgG3	Ig fraction - Protein A
YII-Y042-EX	Anti PACAP (1-15)	50 UL	Human Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y050-EX	Anti PACAP 27	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
YII-Y051-EX	Anti PACAP 27 (13-27)	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
YII-Y040-EX	Anti PACAP 38	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y041-EX	Anti PACAP 38 (22-38)	50 UL	Human Rat	Rabbit		IHC	polyclonal		Serum
SML-ROI002-EX	Anti PAD2	0.1 ML	Human	Rabbit		WB IHC(p)	polyclonal	IgG	Purified
YII-Y090-EX	Anti Pancreastatin (33-51)	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
KAL-KR059	Anti PDX1	25 UG	Mouse	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KH012	Anti Pentosidine	50 UG	-	Mouse	PEN-12	IHC ELISA	monoclonal	IgG1	Affinity Purified
KAL-KH012-01	Anti Pentosidine Biotin	50 UG	-	Mouse	PEN-12	IHC ELISA	monoclonal	IgG1	Affinity Purified
KAL-KH012-02	Anti Pentosidine HRP	50 UG	-	Mouse	PEN-12	IHC ELISA	monoclonal	IgG1	Affinity Purified
KAL-KI044	Anti Period1 (PER1)	200 UG	Human	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
KAL-KI047	Anti Period1 (PER1)	200 UG	Rat	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
KAL-KI046	Anti Period3 (PER3)	200 UG	Human	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
KAL-KI048	Anti Period3 (PER3)	200 UG	Rat	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
ATA-CB-P1RA	Anti Peroxiredoxin 1	0.1 ML	Human Rat	Rabbit		IHC(f) IHC(p) WB ELISA IF LM	polyclonal		Serum
ATA-CB-P2RA	Anti Peroxiredoxin 2	0.1 ML	Human Rat	Rabbit		IHC(f) IHC(p) WB ELISA IF LM	polyclonal		Serum
ATA-CB-P4RA	Anti Peroxiredoxin 4	0.1 ML	Human Rat	Rabbit		IHC(f) IHC(p) WB ELISA IF LM	polyclonal		Serum
ATA-CB-P6RA	Anti Peroxiredoxin 6	0.1 ML	Human Rat	Rabbit		IHC(f) IHC(p) WB ELISA IF LM	polyclonal		Serum
ATA-CB-P6SA	Anti Peroxiredoxin 6	0.1 ML	Human Rat Mouse	Sheep		IHC(f) IHC(p) WB ELISA IF LM	polyclonal		Serum
YII-Y020-EX	Anti PHI	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
YII-Y022-EX	Anti PHI	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y021-EX	Anti PHI (13-27)	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
KAL-KR076	Anti phospho-Tau (Ser416)	25 UG	Human Mouse	Rabbit		IHC	polyclonal	IgG	Affinity Purified
ALF-274042041	Anti PIP2	0.20 MG		Mouse	KT10	ELISA IHC IB WB	monoclonal	IgG2b	Ig fraction - Protein A
KAL-KW207	Anti PLGF	100 UG	Human Mouse Rat	Rabbit		WB IHC(p) IC	polyclonal	IgG	Affinity Purified
PPX-PP-H7223-00	Anti PNR	0.1 ML	Human Rat	Mouse	H7223	WB ELISA IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
KAL-KR052	Anti POEM (Nephronectin)	50 UG	Human	Rabbit		IHC	polyclonal	IgG	Affinity Purified
YII-Y080-EX	Anti PP (18-36)	50 UL	Human Rat	Rabbit		IHC RIA	polyclonal		Serum
PPX-PP-A3409A-00	Anti PPAR $\gamma$ common	0.1 ML	Human Mouse Rat	Mouse	A3409A	WB ELISA IP Gel Shift IHC ChIP	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-H5344-00	Anti PR (isoform B)	0.1 ML	Human Rat Mouse	Mouse	H5344	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-A9621A-00	Anti PR common	0.1 ML	Human Rat Mouse	Mouse	A9621A	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
KAL-KG406	Anti Prei4	25 UG	Mouse	Rabbit		IHC(p)	polyclonal		Affinity Purified
YII-Y185-EX	Anti Preprogalanin (89-124)	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y190-EX	Anti Preprorenin (21-64)	50 UL	Human	Rabbit		IHC(f) RIA	polyclonal		Serum
YII-Y191-EX	Anti Preprorenin (288-317)	50 UL	Human	Rabbit		IHC(f) EIA	polyclonal		Serum
LSL-LB-3117	Anti Prion Protein A	100 UL	Human Mouse Rat Bovine Sheep	Rabbit		WB IHC(p) IHC(f) IF	polyclonal		Serum
LSL-LB-3227	Anti Prion Protein B	100 UL	Human Mouse Rat Bovine Sheep	Rabbit		WB IHC(p) IHC(f) IF	polyclonal		Serum
YII-Y310-EX	Anti ProLH RH (1-26)	50 UL	Human Rat Porcine	Rabbit		IHC(f) RIA	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
YII-Y421-EX	Anti PTH (1-15)	50 UL	Mouse Rat	Rabbit		IHC(p) IHC(f) RIA	polyclonal		Serum
YII-Y420-EX	Anti PTH (1-34) NH2	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y201-EX	Anti PTHrP (1-34)-NH <sub>2</sub>	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
YII-Y202-EX	Anti PTHrP (15-34)-NH <sub>2</sub>	50 UL	Human Rat	Rabbit		IHC(p) IHC(f) RIA	polyclonal		Serum
PPX-PP-PPZ1228-00	Anti PTX3	0.1 ML	Human	Mouse	PPZ1228	WB IHC ELISA	monoclonal	IgG2b	Ig fraction - Ammonium Sulphate
KAL-KH010	Anti Pyrraline	20 UG	-	Mouse	H12	IHC ELISA	monoclonal	IgG1	Affinity Purified
KAL-KH010-01	Anti Pyrraline Biotin	20 UG	-	Mouse	H12	IHC ELISA	monoclonal	IgG1	Affinity Purified
KAL-KH010-02	Anti Pyrraline HRP	20 UG	-	Mouse	H12	IHC	monoclonal	IgG1	Affinity Purified
YII-Y070-EX	Anti PYY	50 UL	Porcine	Rabbit		IHC RIA	polyclonal		Serum
YII-Y072-EX	Anti PYY	50 UL	Human	Rabbit		IHC(p)	polyclonal		Serum
KAL-KH039	Anti RAGE	100 UG	Human	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
YII-Y441-EX	Anti Rattin 24	50 UL		Rabbit		IHC EIA	polyclonal		Serum
YII-Y280-EX	Anti Rb (901-928)	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum
SIM-2ZRBP1	Anti RBP-Jk	0.5 ML	Mouse	Rat	K0043	IHC ELISA IP	monoclonal	IgG1	Purified
SIM-2ZRBP2	Anti RBP-Jk	2 ML	Mouse	Rat	T6709	WB IHC ELISA not_IP	monoclonal	IgG2a	Supernatant
KAL-KE030	Anti Related To b(0,+)-Type Amino Acid Transporter (rBAT)	25 UG	Rat	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
PPX-PP-H2729-00	Anti REVERB $\beta$	0.1 ML	Human Rat	Mouse	H2729	WB ELISA IP IHC	monoclonal	IgG3	Ig fraction - Ammonium Sulphate
LSL-LB-5509	Anti Rhodopsin	100 UL	Frog Squid	Rabbit		IHC(p) IHC(f)	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
LSL-LB-5555	Anti Rhodopsin	100 UL	Human Mouse Rat Bovine Chicken Frog Vertebrate	Rabbit		WB IHC(p) IHC(f)	polyclonal		Serum
PPX-PP-H3925-00	Anti ROR common	0.1 ML	Human Rat Mouse	Mouse	H3925	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
PPX-PP-K8508-00	Anti RXR $\alpha$	0.1 ML	Human Mouse Rat	Mouse	K8508	WB ELISA IP Gel Shift IHC ChIP	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
LSL-LB-9197	Anti S 100 Protein	100 UL	Human Mouse Rat Bovine Goat Eel	Rabbit		IHC(p) IHC(f)	polyclonal		Serum
YII-YP080-EX	Anti S100 $\beta$ (41-60)	50 UL	Human Mouse Rat	Rabbit		WB IHC	polyclonal		Serum
YII-YP081-EX	Anti S100 $\beta$ (74-92)	50 UL	Human Mouse Rat	Rabbit		WB IHC	polyclonal		Serum
KAL-KY008	Anti S19 Ribosomal Protein	100 UG	Human	Rabbit		IHC ELISA	polyclonal	IgG	Affinity Purified
PPX-PP-K9814-00	Anti SALL1	0.1 ML	Human Mouse	Mouse	K9814	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
BAM-65-101	Anti SARS	50 UG		Mouse	3A2	WB IHC ELISA	monoclonal		
BAM-65-102	Anti SARS	250 UG		Mouse	3A2	WB IHC ELISA	monoclonal		
CAC-GU01-M01AS-A	Anti Schwann Cell/peripheral myelin	50 UL	Human Mouse Rat	Mouse	Schwann/2E	IHC	monoclonal	IgM	Ascites
YII-Y030-EX	Anti Secretin	50 UL	Human Rat Monkey	Rabbit		IHC RIA	polyclonal		Serum
YII-Y031-EX	Anti Secretin	50 UL	Human Rat Porcine	Rabbit		IHC RIA	polyclonal		Serum
YII-Y032-EX	Anti Secretin	50 UL	Rat	Rabbit		IHC RIA	polyclonal		Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
KAL-KR100	Anti Semaphorin 3B	100 UG	Human	Rabbit		IHC	polyclonal	IgG	Affinity Purified
CAC-SDT-02-SP5	Anti Septin5, phospho Ser17	100 UL	Mouse	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
PPX-PP-N1665-00	Anti SF-1	0.1 ML	Human Rat Mouse	Mouse	N1665	WB ELISA IP IHC	monoclonal	IgG1	Ig fraction - Ammonium Sulphate
CAC-TNL-002-SH1	Anti SHP-1	100 UL		Rabbit		WB IHC IP	polyclonal	IgG	Serum
CAC-TNL-002-SH2	Anti SHP-1, phospho Tyr536	100 UL		Rabbit		WB IHC IP	polyclonal	IgG	Serum
BAM-73-001	Anti SLA (clone TRA54)	50 UG		Mouse	TRA54	WB IHC			
BAM-73-002	Anti SLA (clone TRA54)	250 UG		Mouse	TRA54	WB IHC			
KAL-KG404	Anti Slc25a19	25 UG	Mouse	Rabbit		IHC(p)	polyclonal		Affinity Purified
YMS-80058	Anti Smooth Muscle Myosin Heavy Chain (SM1)	100 UG	Mouse Human	Rat	KM3669	WB IHC	monoclonal	IgG2a	Purified
SML-ROI001-EX	Anti SMP30 [Regucalcin, Gluconolactonase (GNL)]	0.1 ML	Human Mouse Rat	Rabbit		WB IHC(p)	polyclonal	IgG	Purified
ATA-CB-SOD1RA	Anti SOD1	0.1 ML	Human	Rabbit		IHC(f) IHC(p) WB ELISA IF LM	polyclonal		Serum
ATA-CB-SOD3RA	Anti SOD3	0.1 ML	Human	Rabbit		IHC(f) IHC(p) WB ELISA IF LM	polyclonal		Serum
YII-YP020-EX	Anti Somatostatin	50 UL	Human Rat	Rabbit		IHC RIA	polyclonal		Serum
ALF-274042010	Anti Sphingosine	0.25 MG		Mouse	NHSPH	IHC IB	monoclonal	IgM	Ig fraction - Ion Exchange /Gel Filtration
ALF-274042034	Anti Sphingosine-1-phosphate [S1P]	0.25 MG		Mouse	NHS1P	IHC IB	monoclonal	IgM	Ig fraction - Ion Exchange /Gel Filtration
CAC-NIG-L1-SHA2	Anti SPT16	100 UL		Rabbit		WB ELISA IHC IP ChIP	polyclonal	IgG	Serum
CAC-NIG-L1-SHA1	Anti SSRP1	100 UL		Rabbit		WB ELISA IHC IP ChIP	polyclonal	IgG	Serum

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
KAL-KW216	Anti STAT1	100 UG	Human Mouse Rabbit Rat	Rabbit		WB IHC(p)	polyclonal	IgG	Affinity Purified
YII-Y150-EX	Anti Substance P	50 UL	Human Rat Canine Guinea Pig	Rabbit		IHC RIA	polyclonal		Serum
KAL-KM108	Anti SulFP1/sulf-1	50 UG	Rat	Rabbit		IHC	polyclonal		Affinity Purified
KAL-KM109	Anti SulFP2/sulf-2	25 UG	Rat	Rabbit		IHC	polyclonal		Affinity Purified
ATA-CB-TAURA	Anti Tau	0.1 ML	Human	Rabbit		IHC(f) IHC(p) IF LM	polyclonal		Serum
ATA-CB-TAUSA	Anti Tau	0.1 ML	Human	Sheep		IHC(f) IHC(p) IF LM	polyclonal		Serum
CAC-TIP-TD-P07	Anti TDP-43 (3-12)	100 UL		Rabbit		WB ELISA IHC	polyclonal	IgG	Serum
CAC-TIP-TD-P09	Anti TDP-43 (405-414)	100 UL		Rabbit		WB ELISA IHC	polyclonal	IgG	Serum
CAC-TIP-PTD-P05	Anti TDP-43, phospho Ser403/404	100 UL		Rabbit		WB ELISA IHC	polyclonal	IgG	Serum
CAC-TIP-PTD-P03	Anti TDP-43, phospho Ser409	100 UL		Rabbit		WB ELISA IHC	polyclonal	IgG	Serum
CAC-TIP-PTD-M01	Anti TDP-43, phospho Ser409/410	50 UL		Mouse	11-9	WB ELISA IHC	monoclonal	IgG1	Supernatant
CAC-TIP-PTD-P01	Anti TDP-43, phospho Ser409/410-1	100 UL		Rabbit		WB ELISA IHC	polyclonal	IgG	Serum
CAC-TIP-PTD-P02	Anti TDP-43, phospho Ser409/410-2	100 UL		Rabbit		WB ELISA IHC	polyclonal	IgG	Serum
CAC-TIP-PTD-P04	Anti TDP-43, phospho Ser410	100 UL		Rabbit		WB ELISA IHC	polyclonal	IgG	Serum
KAL-KR111	Anti Tensin 2	50 UG	Human	Rabbit		WB IHC	polyclonal		Affinity Purified

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
YII-Y240-EX	Anti TGF $\alpha$	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
YII-Y241-EX	Anti TGF $\beta$ 1	50 UL	Human	Rabbit		IHC RIA	polyclonal		Serum
DFK-F-23-EX	Anti TIMP 1	1 ML	Bovine Human Rat	Mouse	7-6C1	IHC(p) WB EIA	monoclonal	IgG1 $\kappa$	Purified
DFK-F-26-EX	Anti TIMP 1	1 ML	Human	Mouse	147-6D11	IHC(p) WB	monoclonal	IgG1 $\kappa$	Purified
DFK-F-70-EX	Anti TIMP 2	1 ML	Human Mouse Rat Guinea Pig Rabbit Bovine	Mouse	67-4H11	IHC(p) IHC(f) WB EIA	monoclonal	IgG1 $\kappa$	Purified
CAC-THU-A-TLR2	Anti TLR2	50 UL		Rabbit		WB IHC ELISA IP FC	polyclonal	IgG	Affinity Purified
CAC-THU-A-TLR9	Anti TLR9	50 UL		Rabbit		WB IHC ELISA IP FC	polyclonal	IgG	Affinity Purified
PPX-PP-H6510-00	Anti TLX	0.1 ML	Human Rat Mouse	Mouse	H6510	WB ELISA IP IHC	monoclonal		
YII-YC032-EX	Anti TNF $\alpha$	50 UL	Rat	Mouse		IHC(p) IHC(f) RIA	polyclonal		Serum
YII-YP040-EX	Anti Total Glucagon	50 UL	Human Rat	Rabbit		IHC RIA	polyclonal		-
PPX-PP-H0107B-00	Anti TR4	0.1 ML	Human Mouse Rat	Mouse	H0107B	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
LSL-LB-1084	Anti Tropomyosin	100 UL	Bovine Chicken	Rabbit		WB IHC	polyclonal		Serum
KAL-KM060	Anti TRPM8 (CMR1)	25 UG	Rat	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KM018	Anti TRPV1 (VR-1)	5 UG	Rat	Rabbit		IHC	polyclonal	IgG	Affinity Purified
KAL-KM019	Anti TRPV2 (VRL-1)	5 UG	Rat	Rabbit		IHC	polyclonal	IgG	Affinity Purified

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
PPX-PP-H2804-00	Anti TR $\alpha$	0.1 ML	Human Mouse Rat	Mouse	H2804	WB ELISA IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
HKD-BC23	Anti ubiquitin	400 UL	Bovine	Mouse	49H9	WB IHC IP	monoclonal	IgG1	
KAL-KE094	Anti URAT1	25 UG	Human	Rabbit		WB IHC	polyclonal	IgG	Affinity Purified
YII-Y361-EX	Anti Urocortin	50 UL	Mouse Rat	Rabbit		IHC(p) IHC(f) RIA	polyclonal		Serum
YII-Y360-EX	Anti Urocortin (3-40)	50 UL	Mouse Rat	Rabbit		IHC RIA	polyclonal		Serum
YII-Y362-EX	Anti Urocortin 1	50 UL	Mouse Rat	Rabbit		IHC ELISA	polyclonal	IgG	Serum
YII-Y363-EX	Anti Urocortin 2	50 UL	Mouse	Rabbit		IHC ELISA	polyclonal	IgG	Serum
YII-Y365-EX	Anti Urocortin 2	50 UL	Rat	Rabbit		IHC EIA	polyclonal		Serum
YII-Y364-EX	Anti Urocortin 3	50 UL	Mouse Rat	Rabbit		IHC ELISA	polyclonal	IgG	Serum
KAL-KG403	Anti Usp11	25 UG	Mouse	Rabbit		IHC(p)	polyclonal		Affinity Purified
PPX-PP-H4537-00	Anti VDR	0.1 ML	Human Mouse Rat	Mouse	H4537	WB ELISA IP IHC	monoclonal	IgG2a	Ig fraction - Ammonium Sulphate
YII-Y010-EX	Anti VIP	50 UL	Human Porcine	Rabbit		IHC RIA	polyclonal		Serum
LSL-LB-2007	Anti Vitronectin	100 UL	Bovine	Rabbit		WB IHC	polyclonal		Serum
LSL-LB-2096	Anti Vitronectin	100 UL	Human Mouse Rat Bovine Porcine Canine Sheep Goat	Rabbit		WB IHC	polyclonal		Serum
KAL-KG407	Anti Zmat4	25 UG	Mouse	Rabbit		IHC(p)	polyclonal		Affinity Purified

Catalog Number	Item Name	Quantity	Species Reactivity	Host Species	Clone Number	Applications	Type	Isotype	Form/Purification
YII-Y300-EX	Anti $\alpha$ Endorphin	50 UL	Rat	Rabbit		IHC(f) RIA	polyclonal		Serum
ATA-CB-AS1S	Anti $\alpha$ synuclein	0.1 ML	Human	Sheep		IHC(f) IHC(p) WB IF LM	polyclonal		Serum
ATA-CB-AS2S	Anti $\alpha$ synuclein	0.1 ML	Human	Sheep		IHC(f) IHC(p) WB IF LM	polyclonal		Serum
ATA-CB-AS3S	Anti $\alpha$ synuclein	0.1 ML	Human	Sheep		IHC(f) IHC(p) WB IF LM	polyclonal		Serum
ATA-CB-BCRS1R	Anti $\alpha/\beta$ crystallin	0.1 ML	Bovine Human	Rabbit		IHC(f) IHC(p) IF LM	polyclonal		Serum
ATA-CB-BS1S	Anti $\beta$ synuclein	0.1 ML	Human	Sheep		IHC(f) IHC(p) IF LM	polyclonal		Serum
SID-IDL031	Anti $\beta$ 2 Microglobulin	100 UG	Human	Mouse	EMR B12	WB IHC	monoclonal	IgG1	
SID-IDL034	Anti $\beta$ 2 Microglobulin Biotin	50 UG	Human	Mouse	EMR B12	WB IHC	monoclonal	IgG1	
HKD-AB47	Anti $\beta$ 2-Microglobulin	100 UL	Human	Mouse	EMRB6-12	WB IHC	monoclonal	IgG1	
ACB-GGT-AGT1EX	Anti $\gamma$ Glutamyl Transpeptidase	100 UG	Human Mouse Rat	Mouse	AGT1	WB IHC	monoclonal	IgG1	Ig fraction - Protein A
SML-ROA002-EX	Anti $\gamma$ subnit, Fc reseptor	200 UL	Mouse Rat	Rabbit		WB IHC(p) IP	polyclonal	IgG	Affinity Purified
ATA-CB-GS1S	Anti $\gamma$ synuclein	0.1 ML	Human	Sheep		IHC(f) WB	polyclonal		Serum