

POLYCLONAL ANTIBODY

For research use only, Not for diagnostic use.

Catalog No. EXO-AB-01

Anti Bovine Milk Exosome

Product type	Primary antibody
Immunogen	Bovine
Raised in	Rabbit
Myeloma	-
Clone number	-
Isotype	-
Source	Serum
Purification	Ammonium sulfate precipitation Milk Whey protein absorbed
Buffer	PBS(-) 15ppm Proclin300
Concentration	1 mg / mL
Volume	100 ul
Label	Unlabeled
Specificity	Whole Exosome of bovine Milk
Cross reactivity	Bovine
Storage	Store below - 20°C Aliquot to avoid cycles of freeze/thaw.

Application notes • Western blotting: 1/1000

Recommended
dilutions

References -

RELATED PRODUCT:

Product Name	Quantity	Maker	Cat#
Milk exosome ,Bovine	1 SET (100uLx10 vial)	CSR	EXBM100L
Milk exosome ,Bovine	1 SET (1mLx10 vial)	CSR	EXBM1000L
Milk Exosome ELISA Kit ,Bovine	1 KIT (96test)	CSR	EXBMEL



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

www.cosmobio.co.jp