

MONOCLONAL ANTIBODY

For research use only, Not for diagnostic use.

Catalog No. KYU-HT-M001

Anti Morphine

BACKGROUND

Morphine is a alkaloid isolated from Papaver somniferum, which is psychoactive.

Product type	Primary antibody
Immunogen	Morphine-BSA
Raised in	Mouse (Balb/c)
Myeloma	SP2/0
Clone number	IB1
Isotype	IgG1 λ
Source	Serum free culture supernatant
Purification	Protein G
Buffer	PBS containing 0.02% NaN_3 as a preservative
Concentration	1mg / mL
Volume	100 UL
Label	Unlabeled
Specificity	morphine
Cross reactivity	Specific to morphine
Storage	Store below -20 degrees (below -70 degrees for prolonged storage) Aliquot to avoid cycles of freeze / thaw.

Application notes • **ELISA** (Ref.1)

Recommended dilutions
Other applications have not been tested.
Optimal dilutions/concentrations should be determined by the end user.

References
1) Shoyama et al., DIRECT DETERMINATION OF OPIUM ALKALOID-BOVINE SERUM-ALBUMIN CONJUGATE BY MATRIX-ASSISTED LASER-DESORPTION IONIZATION MASS-SPECTROMETRY. Biol Pharm Bull. 1993 Oct;16(10):1051-3.
PMID: [8287039](#)

For research use only, Not for diagnostic use.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619