

MONOCLONAL ANTIBODY

For research use only, Not for diagnostic use.

Catalog No. KYU-HT-M002

Anti Thebaine

BACKGROUND

Thebaine is a alkaloid isolated from Papaver somniferum.

Product type	Primary antibody
Immunogen	Thebaine-BSA
Raised in	Mouse (Balb/c)
Myeloma	SP2/0
Clone number	5F7
Isotype	IgG2b λ
Source	Serum free culture supernatant
Purification	Protein G
Buffer	PBS containing 0.02% NaN ₃ as a preservative
Concentration	1mg / mL
Volume	100 UL
Label	Unlabeled
Specificity	thebaine
Cross reactivity	Specific to thebaine
Storage	Store below -20 degrees (below -70 degrees for prolonged storage) Aliquot to avoid cycles of freeze / thaw.

Application notes • **ELISA** (Ref.1)

Recommended dilutions
Other applications have not been tested.
Optimal dilutions/concentrations should be determined by the end user.

References
1) Shoyama et al., DIRECT DETERMINATION OF OPIUM ALKALOID-BOVINE SERUM-ALBUMIN CONJUGATE BY MATRIX-ASSISTED LASER-DESORPTION IONIZATION MASS-SPECTROMETRY. Biol Pharm Bull. 1993 Oct;16(10):1051-3. PMID: [8287039](https://pubmed.ncbi.nlm.nih.gov/8287039/)

For research use only, Not for diagnostic use.



COSMO BIO Co., LTD.

[JAPAN]
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

[Outside Japan]
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com