

MONOCLONAL ANTIBODY

For research use only, Not for diagnostic use.

Catalog No. KYU-HT-M003

Anti Ginsenoside Rb1

BACKGROUND

Ginsenoside Rb1 is a saponin isolated from Panax ginseng.

Product type	Primary antibody
Immunogen	Ginsenoside Rb1-BSA
Rased in	Mouse (Balb/c)
Myeloma	SP2/0
Clone number	9G7
Isotype	IgG κ
Source	Serum free culture supernatant
Purification	Protein G
Buffer	PBS containing 0.02% NaN ₃ as a preservative
Concentration	1mg / mL
Volume	100 UL
Label	Unlabeled
Specificity	Ginsenoside Rb1
Cross reactivity	-
Storage	Store below -20 degrees (below -70 degrees for prolonged storage) Aliquot to avoid cycles of freeze / thaw.

Application notes • **ELISA** (Ref.1)

Recommended • **WB** (Ref.2)

dilutions

Other applications have not been tested.

Optimal dilutions/concentrations should be determined by the end user.

References

- 1) Tanaka H. et al., Formation of monoclonal antibody against a major ginseng component, ginsenoside Rb1 and its characterization, Cytotechnology, 29,115, 1999. PMID:[22359060](#)
- 2) Fukuda N, Tanaka H, Shoyama Y., Double staining of ginsenosides by Western blotting using anti-ginsenoside Rb1 and Rg1 monoclonal antibodies. Biol Pharm Bull. 2001 Oct;24(10):1157-60. PMID:[11642323](#)

For research use only, Not for diagnostic use.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619