

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. KZ-HS-P13

Anti Progesterone

BACKGROUND

Progesterone is produced in the ovaries (to be specific, after ovulation in the corpus luteum), the adrenal glands (near the kidney), and, during pregnancy, in the placenta. Progesterone is also stored in adipose (fat) tissue. In humans, increasing amounts of progesterone are produced during pregnancy.

Product type Primary antibodies

HostRabbitFormSerumVolume50 μl

Antigen Progesterone 3-CMO:BSA

Cross reactivity Human, Rat, Mouse, Bovine, Chicken

Application notes Recommended dilutions

RIA: 1/100,000 - 1/500,000

EIA: n/d

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References 1) 長谷川喜久 (1988) 内分泌試験 (各種ホルモンの抗体作製法) 佐久間勇次 監修

ウサギー生殖生理と実験手技 P.201-216、近代出版

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.

[JAPAN]

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU. TOKYO 135-0016, JAPAN Phone: +81-3-5632-9610

Phone: +81-3-5632-9610 FAX: +81-3-5632-9619

URL: https://www.cosmobio.co.jp/



COSMO BIO USA

[Outside Japan] 2792 Loker Ave West, Suite 101 Carlsbad, CA 92010, USA email: info@cosmobiousa.com Phone/FAX: (+1) 760-431-4600 URL: www.cosmobiousa.com