

MONOCLONAL ANTIBODY

For research use only, Not for diagnostic use.

Catalog No.NU-01-KE8

Anti keratin 8 chain monoclonal antibody

Product type Primary Antibodies

Immunogen TritonX-100 insoluble cytoskeletal fraction of rabbit intestinal epithelial cells.

Clone number RL273

Isotype IgG1, kappa

Host Mouse

Formulation Hybridoma supernatant with 0.02% NaN3 as a preservative.

Volume 500ul

LabelUnlabeledSpecificityKeratin 8Cross reactivityHuman, rabbit

Storage Store at -20° C or -70° C in small aliquots for prolonged storage.

Repeated freeze-thaw cycles can damage the immunoreactivity of an antibody.

Application notes

Recommended use WB, IF

Not tested yet in other applications.

Recommended Western Blot: 1:100–1:500 for detection of a 50-kDa polypeptide in whole cell lysate of DLD-1 cells.

dilutions Immunofluorescence: 1:100–1:250 for staining of acetone-fixed cryostat frozen tissue sections.

Optimal dilutions must be determined by end user.




Fig.1 Western blot analysis

Whole cell lysate prepared from DLD-1 cells was immunoblotted with RL273 antibody (1:200 dilution).

The antibody detected a polypeptide band for keratin 8 at approximately 50 kDa.

DLD-1 is a human colorectal adenocarcinoma cell line.




Fig.2 Immunofluorescence microscopy of rabbit colon and DLD-1 cells.

A frozen section of rabbit colon was stained by RL273 antibody at 1:250 dilution (A). The absorptive epithelial cells and mucosal gland cells were decorated by the antibody. Panel B shows the phase contrast image of A. Bar = 100um. Acetone-fixed DLD-1 cells were stained by RL273 antibody at 1:250 dilution (C), showing the filamentous networks of keratins (green). Cell nuclei were visualized by DAPI (blue). Bar = 25um.

