



POLYCLONAL ANTIBODY

Catalog No. SDT-01-AO1

Anti-AAO1 (Arabidopsis Aldehyde Oxidase 1)

Product type	Primary antibodies
Host	Rabbit
Source	Serum (Partially purified (0-40% ammonium sulfate fraction))
Form	TBS (pH8.0) with 0.02% NaN ₃ as a preservative.
Volume	200ul
Concentration	
Specificity	Cross reacts with other AAOs in WB. Not cross reacts in IP.
Antigen	Arabidopsis thaliana aldehyde oxidase 1 (peptide corresponding to 288 amino acids (1081-1368aa))
Isotype	IgG

Application notes WB, IP IHC, ELISA, FC: not tested

Recommended use

Recommended dilutions

Western Blot: 1/1,000. Predicted molecular weight: 150 kDa

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Plant (Arabidopsis, Pea etc.)

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References

Akaba, S., Seo, M., Dohmae, N., Takio, K., Sekimoto, H., Kamiya, Y., Furuya, N., Komano, T. and Koshiba, T. (1999) Production of homo- and hetero-dimeric isozymes from two aldehyde oxidase genes of Arabidopsis thaliana. J Biochem (Tokyo) 126:395-401.

This product was made by using the research results of Tokyo Metropolitan University.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

[JAPAN]
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

[Outside Japan]
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com