



COSMO BIO CO., LTD.

Inspiration for Life Science

POLYCLONAL ANTIBODY

Catalog No. SDT-01-AO3

Anti-AAO3 (Arabidopsis Aldehyde Oxidase 3)

Product type	Primary antibodies
Host	Rabbit
Source	Serum (Partially purified (0-40% ammonium sulfate fraction))
Form	TBS (pH8.0) with 0.02% NaN ₃ as a preservative.
Volume	200ul
Concentration	
Specificity	Cross reacts with other AAOs in WB. Not cross reacts in IP.
Antigen	Arabidopsis thaliana aldehyde oxidase 3 (peptide corresponding to 411 amino acids (347-757aa))
Isotype	IgG

Application notes WB, IP IHC, ELISA, FC: not tested

Recommended use

Recommended dilutions

Western Blot: 1/1,000. Predicted molecular weight: 130 kDa

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Plant (Arabidopsis, Pea etc.)

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References Seo, M., Koiwai, H., Akaba, S., Komano, T., Oritani, T., Kamiya, Y. and Koshiba, T. (2000) Absciscic aldehyde oxidase in leaves of Arabidopsis thaliana. Plant J. 23:481-488.

This product was made by using the research results of Tokyo Metropolitan University.

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com