



POLYCLONAL ANTIBODY

Catalog No. SDT-01-APX

Anti-Cytosolic Ascorbate Peroxidase

Product type	Primary antibodies
Host	Rabbit
Source	Serum (Partially purified (0-40% ammonium sulfate fraction))
Form	TBS (pH7.5) with 0.02% NaN ₃ as a preservative.
Volume	200ul
Concentration	
Specificity	Cytosolic Ascorbate Peroxidase
Antigen	Cytosolic ascorbate peroxidase purified from Zea mays
Isotype	IgG

Application notes WB, IHC IP, ELISA, FC: not tested

Recommended use

Recommended dilutions

Western Blot: 1/1,000-1/2,000. Predicted molecular weight: 145 kDa 28kDa

Immunohistochemistry: 1/100-1/500

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Plant (Maize, Rice, etc.)

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References Koshiha, T. (1993) Cytosolic ascorbate peroxidase in seedlings and leaves of maize (Zea mays). Plant Cell Physiol. 34:713-721.

This product was made by using the research results of Tokyo Metropolitan University.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

【JAPAN】
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com