



POLYCLONAL ANTIBODY

For research use only, Not for diagnostic use.

Catalog No. SDT-02-P35

Anti p35, phospho Ser 8

Product type	Primary Antibodies
Host	Rabbit
Purification	Immunogen affinity purified
Form	PBS (pH7.4) with 0.1% NaN ₃ as a preservative.
Volume	100 µl
Immunogen	Phosphorylated Ser8 of p35 Cdk5 activator
Cross reactivity	Human, Mouse, Rat
Storage	Store below -20°C (below -70°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.

Application notes	• Western Blot: 1/200. Predicted molecular weight: 35kDa
Recommended dilutions	• IHC, IP, ELISA, FC: not tested. 1) Kamei H., <i>et al.</i> , J. Biol. Chem. 282(3): 1687-1694 (2007), PMID: 17121855
References	

This product was made by using the research results of Tokyo Metropolitan University.

For research use only, Not for diagnostic use.



COSMO BIO Co., LTD.

【JAPAN】
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com