



Anti p35, phospho Ser 8

Catalog No. SDT-02-P35

Lot number	*****
Host animal	Rabbit
Immunogen	Phosphorylated Ser8 of p35 Cdk5 activator
Purity	Immunogen affinity purified
Volume	100 µl
Formulation	PBS (pH7.4) with 0.1% NaN ₃ as a preservative.
Cross Reactivity	Human, Mouse, Rat
Application/ Recommended concentration	<ul style="list-style-type: none">• Western Blot: 1/200. Predicted molecular weight: 35kDa• IHC, IP, ELISA, FC: not tested yet.
Storage	Store below -20°C (below -70°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.
Reference	Kamei, H., Saito, T., Ozawa, M., Fujita, Y., Asada, A., Bibb, J.A., Saido, T.C., Sorimachi, H. and Hisanaga, S. Suppression of Calpain-dependent Cleavage of the CDK5 Activator p35 to p25 by Site-specific Phosphorylation <i>J. Biol. Chem.</i> 282(3): 1687-1694 (2007)

For research use only; not for use as a diagnostic.

This product was made by using the research results of Tokyo Metropolitan University.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[International/Outside Japan] Phone: +81-3-5632-9617

[国内連絡先] Phone: 03-5632-9610

FAX: +81-3-5632-9618

FAX: 03-5632-9619