

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. SK-T01-002

Anti-PACE4 propeptide antiserum

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid
Volume	100 µl
Concentration	
Specificity	
Antigen	Recombinant human PACE4 (SPC4) propeptide (amino acid number 69-149)-DHFR fusion protein
Isotype	IgG

Application notes WB, IHC, IP, ELISA
Recommended use

Recommended dilutions

Western blotting: 1/5,000

Predicted molecular weight: Precursor 110 kDa, Mature form 105 kDa

Immunohistochemistry: 1/5,000 to 1/10,000

Immunoprecipitation/ELISA: 1/1,000

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Not tested

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References 1) Taniguchi et al (2002) A critical role for the carboxy terminal region of the proprotein convertase, PACE4A, in the regulation of its autocatalytic activation coupled with secretion. Biochem. Biophys. Res. Commun. 290, 878-884

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619