

MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. SU-IZ-M01

Anti human SCCA2

BACKGROUND

Anti-SCCA2 antibody is a rat monoclonal antibody which obtained from the immunization with purified, E. coli-derived, recombinant human SCCA2. This antibody can be used for the detection of human SCCA2 by immunoprecipitation and ELISA with no cross-reaction to SCCA1.

Product type Primary antibodies

Host Rat

Source Culture supernatant

Form Liquid

Purified IgG with 0.05% NaN_3 as a preservative.

Volume 100 ul
Concentration 0.54 mg/,l
Specificity SCCA2

Antigen Recombinant SCCA2 expressed by E. coli.

Clone SS8J Isotype IgG2b/k

Application notes IP, ELISA

Recommended dilutions

IP, 1-10 ug/ml ELISA, 1-10 ug/ml

Optimal dilutions/concentrations should be determined by the end user.

Cross reactivity Human

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp e-mail: export@cosmobio.co.jp

[Outside Japan] Phone: +81-3-5632-9617 [国内連絡先] Phone: +81-3-5632-9610 FAX: +81-3-5632-9618 FAX: +81-3-5632-9619