

**MONOCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. SU-IZ-M03**

## Anti human IL-13Ra D1 domain (SS4B)

### BACKGROUND

Anti-SCCA2 antibody is a rat monoclonal antibody which obtained from the immunization with purified, E. coli-derived, recombinant human SCCA2. This antibody can be used for the detection of human SCCA2 by immunoprecipitation and ELISA with no cross-reaction to SCCA1.

<b>Product type</b>	Primary antibodies
<b>Host</b>	Rat
<b>Source</b>	Mouse ascites
<b>Form</b>	Liquid
	Purified IgG with 0.05% NaN <sub>3</sub> as a preservative.
<b>Volume</b>	100 ul
<b>Concentration</b>	0.47 mg/ml
<b>Specificity</b>	IL-13 receptor $\alpha$ 1D1 domain
<b>Antigen</b>	Recombinant human D1 domain expressed by E. coli
<b>Clone</b>	SS4B
<b>Isotype</b>	IgG2a/k

<b>Application notes</b>	FC
	<b>Recommended dilutions</b>
	FC, 1-10 ug/ml
	Optimal dilutions/concentrations should be determined by the end user.
<b>Cross reactivity</b>	Human
<b>Storage</b>	Store below -20°C (below -70°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.

*For research use only. Not for clinical diagnosis.*



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619