

MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. SU-IZ-M07

Anti human SCCA1/SCCA2

BACKGROUND

Anti-SCCA1/2 antibody is a rat monoclonal antibody which obtained from the immunization with purified, E. coli-derived, recombinant human SCCA1. This antibody can be used for the detection of human SCCA1 and SCCA2 by immunoprecipitation and ELISA.

Product type Primary antibodies

Host Rat

Source Culture supernatant

Form Liquid

Purified IgG with 0.05% NaN₃ as a preservative.

Volume 100 ul Concentration 0.3 mg/,l

Specificity SCCA1, SCCA2

Antigen Recombinant SCCA1 expressed by E. coli.

Clone SS6A Isotype IgG2b/k

Application notes IP, ELISA

Recommended dilutions

IP, 1-10 ug/ml ELISA, 1-10 ug/ml

Optimal dilutions/concentrations should be determined by the end user.

Cross reactivity Human

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp e-mail: export@cosmobio.co.jp

[Outside Japan] Phone: +81-3-5632-9617 [国内連絡先] Phone: +81-3-5632-9610 FAX: +81-3-5632-9618 FAX: +81-3-5632-9619