

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. SU-IZ-P01

Anti EGFP

BACKGROUND

Anti EGFP antibody is a rabbit polyclonal antibody which obtained from the immunization with purified, E. coli-derived, recombinant EGFP. This antibody has been tested and confirmed for use in immunoblotting against enhanced GFP.

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid
	Serum with 0.05% NaN ₃ as a preservative.
Volume	100 ul
Concentration	
Specificity	GFP
Antigen	Recombinant EGFP expressed by E. coli.

Application notes

WB

Recommended dilutions

Western blotting, 1/1,000 to 1/5,000

Optimal dilutions/concentrations should be determined by the end user.

Storage

Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619