



**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. SU-IZ-P04**

## Anti human SCCA

### BACKGROUND

Anti-SCCA antibody is a rabbit polyclonal antibody which obtained from the immunization with purified, E. coli-derived, recombinant human SCCA2. This antibody can be used for the detection of SCCA by immunoblotting and immunostaining.

<b>Product type</b>	Primary antibodies
<b>Host</b>	Rabbit
<b>Source</b>	Serum
<b>Form</b>	Liquid
	Serum with 0.05% NaN <sub>3</sub> as a preservative.
<b>Volume</b>	100 ul
<b>Specificity</b>	SCCA1, SCCA2
<b>Antigen</b>	Recombinant SCCA2 expressed by E. coli.

---

<b>Application notes</b>	WB, IHC <b>Recommended dilutions</b> Western blotting, 1/100 to 1/1,000 Immunohistochemistry, 1/50 to 1/500 Optimal dilutions/concentrations should be determined by the end user.
<b>Cross reactivity</b>	Human
<b>Storage</b>	Store below -20°C (below -70°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.
<b>References</b>	1) Mitsuishi, K. et al., Clin. Exp. Allergy. 35:1327-1333, 2005