

## POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

## Catalog No. SU-IZ-P05

## Anti filamin A C-terminal region

## BACKGROUND

Anti-filamin A antibody is a rabbit polyclonal antibody which obtained from the immunization with purified, E. coli-derived, recombinant mouse filamin A. This antibody can be used for the detection of filamin A by immunoblotting.

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid
	Serum with 0.05% $NaN_3$ as a preservative.
Volume	100 ul
Concentration	
Specificity	Filamin A
Antigen	Recombinant mouse filamin A expressed by E. coli.

Application notes	WB
	Recommended dilutions
	Western blotting, 1/500 to 1/2,000
	Optimal dilutions/concentrations should be determined by the end user.
Cross reactivity	Mouse
Storage	Store below -20°C (below -70°C for prolonged storage).
	Aliquot to avoid cycles of freeze/thaw.

For research use only. Not for clinical diagnosis.

