



POLYCLONAL ANTIBODY

Catalog No. YCU-PS-M2

Anti muscarinic acetylcholine receptor m2

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid
Volume	1 ml
Concentration	
Specificity	Muscarinic acetylcholine receptor m2
Antigen	Muscarinic acetylcholine receptor m2 (N227-C324)
Isotype	

Application notes IP Not yet tested in other applications.

Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Human, Porcine

Storage Store below -20°C (below -70°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.

References 1) K. Shiozaki, E. Iseki, H. Uchiyama, Y. Watanabe, T. Haga, K. Kameyama, T. Ikeda, T. Yamamoto, K. Kosaka Alterations of muscarinic acetylcholine receptor subtypes in diffuse Lewy body disease: relation to Alzheimer's disease J Neurol Neurosurg Psychiatry. 1999 Aug;67(2):209-13.

2) Kazumasa Shiozaki, Eizo Iseki, Hiroaki Hino, Kenji Kosaka Distribution of m1 muscarinic acetylcholine receptors in the hippocampus of patients with Alzheimer's disease and dementia with Lewy bodies-an immunohistochemical study J Neurol Sci. 2001 Dec 15;193(1):23-8.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

【JAPAN】
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU, TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com