



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. COP-080055

Anti glyoxysomal lipase

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid
Volume	100 µl
Concentration	Less than 0.1% NaN ₃ as a preservative.
Specificity	Glyoxysomal lipase
Antigen	Purified glyoxysomal lipase
Isotype	

Application notes

WB

Recommended use

Recommended dilutions

Western Blot: ~1/2,000

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Glyoxysomal lipase (castor bean, *Ricinus communis* L.) (62 kDa)
Also known as alkaline lipase

Storage

Store below -20°C
Aliquot to avoid cycles of freeze/thaw.

References

- 1) Maeshima, M. and Beevers, H. (1985) Purification and properties of glyoxysomal lipase from castor bean. *Plant Physiology*, 79: 489-493.

Manufactured by Operon Biotechnologies, K.K



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619