

## **POLYCLONAL ANTIBODY**

For research use only. Not for clinical diagnosis.

Catalog No. COP-080062

## Anti Calpain (CAPN1, 2, 3, 8)

**Product type** Primary antibodies

**Host** Mouse

Source ascitic fluid (expanded with PBS)

 Clone
 1D10A7

 Volume
 50 μl

Immunogen fusion recombinant protein: NP\_001739+NP\_001003962 full length

Speciaficity CAPN1, 2, 3, 8

Cross reactivity Human

Application notes Recommended use: WB

Recommended dilutions
•Western Blot: 1/2,500-1/5,000

Optimal dilutions/concentrations should be determined by the end user.

Store below -20°C (below -80°C for prolonged storage).

**Storage** Aliquot to avoid cycles of freeze/thaw.

**Reference** Kasai Y et al J Biochem 100, 183-190, 1986

## Detection of CAPN1 and CAPN2 by anti-pan-Calpain (clone 1D10A7)

