



Anti CHICO

Product type	Primary Antibody
Immunogen	fusion recombinant protein : Drosophila CHICO (first Met to the 666th Phe)
Host Species	Rabbit
Clone Designation	-
Isotype	-
Source	-
Purification	-
Formulation	Serum
Concentration	-
Volume	100ul
Label	Unlabeled
Specificity	-
Cross reactivity	-
Storage	Store below -20°C (below -80°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.

Application notes	• Western Blot: 1/2,000
Recommended dilutions	• Immunohistochemistry: 1/400

References	Neurosci Res. 2012 May;73(1):49-55. doi: 10.1016/j.neures.2012.02.001. Epub 2012 Feb 10. Mutations in the Drosophila insulin receptor substrate, CHICO, impair olfactory associative learning. Fig : Fig1C, Fig2A
-------------------	---



COSMO BIO CO., LTD.

2-20, Toyo 2-Chome, Koto-ku, Tokyo 135-0016, JAPAN

T: +81-3-5632-9617

日本のお客様

F: +81-3-5632-9618

T: 03-5632-9610

E: export@cosmobio.co.jp

F: 03-5632-9619

W: www.cosmobio.com

W: www.cosmobio.co.jp