



## Anti-hTIMP-1 monoclonal antibody purified mouse IgG

This is a purified mouse monoclonal antibody to recombinant human tissue inhibitor of metalloproteinases-1 (r-hTIMP-1).

Product No. DFK-F-26-EX  
 Clone No. 147-6D11  
 Subclass / type IgG1 / kappa  
 Specificity / Immunogen The antibody specifically reacts with hTIMP-1.

Species	human	mouse	rat	rabbit	Other
Reactivity	(+)	(-)	(-)	NT	NT

NT; not tested

Applications (reference concentration)	Western Blotting	1 µg / mL
	Paraffin Sections	5 - 10 µg / mL <sup>3)</sup>
	Frozen Sections	NT

NT; not tested Optimal working dilutions must be determined by end user.

Concentration 500 µg / mL  
 Presentation Liquid in 0.1M Na-Phosphate buffer, pH7.0, with 0.5 % protease-free BSA and 0.05 % sodium azide.  
 Volume 1 mL  
 Storage Maintain at 2-8 °C in undiluted aliquots for up to 12 months. Do not freeze.  
 References 1) Kodama, S. et al., Collagen Rel. Res., 7, 341-350, 1987.  
 2) Okada, Y. et al., Ann. Rheum. Dis., 48, 645-653, 1989.  
 3) Umehara, F. et al., J. Neuropathol. Exp. Neurol., 57, 839-849, 1998.

Note For laboratory use only

Sep., 2007



**COSMO BIO CO., LTD.**  
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

Phone : +81-3-5632-9617

e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

FAX : +81-3-5632-9618