



Anti-hMMP-13 monoclonal antibody purified mouse IgG

This is a purified mouse monoclonal antibody to recombinant human matrix metalloproteinase 13 (hMMP-13, procollagenase-3).

Product No. DFK-F-94-EX
Clone No. 181-14G11
Subclass / type IgG1 / kappa
Specificity / The antibody specifically reacts with latent and active forms of hMMP-13, and does not
Immunogen cross-react with hMMP-1, 2, 3, 7, 8, 9.

Species	human	mouse	rat	rabbit	Other
Reactivity	(+)	(-) ⁵⁾	NT	NT	NT

NT; not tested

Applications (reference concentration)	Western Blotting	1 µg / mL
	Paraffin Sections	2.5 µg / mL ⁴⁾
	Frozen Sections	NT

NT; not tested Optimal working dilutions must be determined by end user.

Concentration 500 µg / mL
Presentation Liquid in 0.1M Na-Phosphate buffer, pH7.0, with 0.5 % protease-free BSA and 0.05 % sodium azide.
Volume 1 mL
Storage Maintain at 2-8 °C in undiluted aliquots for up to 12 months. Do not freeze.
References 1) Freije, J. M. P. et al., J. Biol. Chem., 269, 16766-16773, 1994.
2) Knäper, V. et al., J. Biol. Chem., 271, 1544-1550, 1996.
3) Knäper, V. et al., J. Biol. Chem., 271, 17124-17131, 1996.
4) Uria, J. A. et al., Am. J. Pathol., 153, 91-101, 1998.
5) Furuichi, I., Internal report of Fuji Chem. Ind. Ltd., 1999.

Note For laboratory use only

Sep., 2007



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

Phone : +81-3-5632-9617

e-mail : export@cosmobio.co.jp

FAX : +81-3-5632-9618