



GeneStem



COSMO BIO CO., LTD.
Inspiration for Life Science

SPECIFICATION SHEET

Description: Anti PSF1 monoclonal

Catalog #: GS001-100, GS001-50

Immunogen: mouse PSF1 **Host:** Rat

Clone#: aho57.2

Subclass: IgG

Concentration: 1.0mg/ml

Cross Reactivity: PSF1, human 100%

Form: Purified by Protein G / Sep4B, supplied in phosphate buffered saline containing 0.09% NaN₃

Storage: store at -20°C.

Applications: Suitable for use in immunohistochemistry, Immunocytochemistry and Western blotting.

Additional Remarks:

Product without NaN₃ is also available on your request.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

11/06/01



COSMO BIO CO., LTD.

2-20, Toyo 2-Chome, Koto-ku, Tokyo 135-0016, JAPAN
T: +81-3-5632-9617 日本のお客様
F: +81-3-5632-9618 T: 03-5632-9610
E: export@cosmobio.co.jp F: 03-5632-9619
W: www.cosmobio.com W: www.cosmobio.co.jp