

SPECIFICATION SHEET

Description:	Anti SLD5 monoclonal		
Catalog #:	GS005-100, GS005-50		
Immunogen:	mouse SLD5	Host:	Rat
Clone#:	AFS3-3-4		
Subclass:	IgG		
Concentration:	1.0mg/ml		
Cross Reactivity:	SLD5, human 100%		
Form:	Purified by Protein G / Sep4B, supplied in phosphate buffered saline containing 0.09% NaN ₃		
Storage:	store at -20°C.		
Applications:	Suitable for use in immunohistochemistry, Immunocytochemistry and Western blotting.		

Additional Remarks:

Product without NaN₃ is also available on your request.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

11/06/01



COSMO BIO CO., LTD.

2-20, Toyo 2-Chome, Koto-ku, Tokyo 135-0016, JAPAN
T: +81-3-5632-9617 | 日本のお客様
F: +81-3-5632-9618 | T: 03-5632-9610
E: export@cosmobio.co.jp | F: 03-5632-9619
W: www.cosmobio.com | W: www.cosmobio.co.jp