



Specification Sheet

Purified HBs Antigen

Code	6Z11
Manufacturing Process	HBsAg positive human plasma
Purification	Potassium bromide equilibrium density-gradient centrifugation. Sucrose density-gradient rate zonal centrifugation.
Buffer	Saline (containing 0.1% NaN ₃)
Protein Concentration	1 mg/mL
Package	1 mL
Storage	Below -20°C
Notice	This product is not inactivated by heat-treatment at 60°C for 10 hours. It should be handled as if it were capable of transmitting, because it is derived from human plasma.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN
<http://www.cosmobio.com> e-mail : export@cosmobio.co.jp
Phone : +81-3-5632-9617 FAX : +81-3-5632-9618